
Midmark M11 - v1.0.0
Total Cycles: 868

Name: RITTER 1
Sterilizer ID: V2178198
06/07/2021 01:40 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 181.3 F 0.0
2:00 181.6 F 3.2
4:00 214.5 F 5.9
6:00 218.1 F 4.0
8:00 226.4 F 4.6
10:00 231.9 F 7.3
12:00 232.9 F 7.3
14:00 246.9 F 13.6

STERILIZING
mm:ss Degrees PSI
0:00 251.6 F 15.9
1:00 252.1 F 16.0
2:00 252.6 F 16.4
3:00 251.8 F 16.2
4:00 252.3 F 16.2
5:00 251.9 F 16.2
6:00 252.4 F 16.3
7:00 251.6 F 15.9
8:00 252.7 F 16.5
9:00 252.2 F 16.2
10:00 251.6 F 15.9
11:00 252.7 F 16.5
12:00 252.4 F 16.3
13:00 251.9 F 16.0
14:00 251.9 F 16.2
15:00 252.7 F 16.5
16:00 252.5 F 16.3
17:00 252.0 F 16.1
18:00 251.6 F 15.9
19:00 252.7 F 16.6
20:00 252.6 F 16.4
21:00 252.4 F 16.2
22:00 251.9 F 16.0
23:00 251.6 F 15.9
24:00 252.6 F 16.6
25:00 252.7 F 16.5
26:00 252.4 F 16.3
27:00 252.1 F 16.1
28:00 251.6 F 15.9
29:00 252.3 F 16.5
30:00 252.7 F 16.5

Min 251.6 F 15.8
Max 252.8 F 16.6

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 199.4 F 0.0
5:00 217.6 F 0.0
10:00 216.3 F 0.0
15:00 219.6 F 0.0
20:00 226.9 F 0.0
25:00 215.6 F 0.0
Dry Cycle Aborted

FILLING: 1:39
HEATING: 15:16
STERILIZING: 30:00
VENTING: 2:56

DRYING: 26:17
TOTAL CYCLE: 01:16:09

06/07/2021 02:56 PM
ACCEPTED: X
REJECTED:
Operator: U01