

-----  
Midmark M11 - v1.0.5  
Total Cycles: 227

Name: RITTER2  
Sterilizer ID: V2313228  
07/20/2021 10:32 AM

BEGIN PACKS CYCLE  
Temp: 250 Degrees F  
Time: 30 Minutes  
Dry: 30 Minutes

FILLING CHAMBER

HEATING  
mm:ss Degrees PSI  
0:00 111.6 F 0.0  
2:00 112.6 F 0.8  
4:00 123.0 F 2.9  
6:00 150.8 F 6.0  
8:00 190.4 F 1.0  
10:00 204.1 F 4.4  
12:00 216.9 F 5.3  
14:00 223.2 F 5.6  
16:00 224.6 F 4.3  
18:00 235.4 F 9.1  
20:00 245.6 F 14.1  
22:00 252.7 F 17.9

STERILIZING  
mm:ss Degrees PSI  
0:00 251.8 F 17.5  
1:00 251.8 F 17.6  
2:00 251.7 F 17.6  
3:00 251.8 F 17.5  
4:00 252.1 F 17.6  
5:00 252.2 F 17.8  
6:00 251.7 F 17.8  
7:00 252.0 F 17.6  
8:00 252.2 F 17.9  
9:00 251.9 F 17.5  
10:00 252.3 F 17.9  
11:00 251.8 F 17.5  
12:00 252.2 F 17.9  
13:00 252.2 F 17.7  
14:00 251.8 F 17.8  
15:00 252.3 F 17.8  
16:00 251.8 F 17.5  
17:00 252.4 F 17.9  
18:00 252.3 F 17.8  
19:00 251.7 F 17.5  
20:00 252.3 F 17.9  
21:00 252.3 F 17.8  
22:00 251.7 F 17.5  
23:00 252.0 F 17.9  
24:00 252.3 F 17.8  
25:00 251.9 F 17.6  
26:00 251.6 F 17.8  
27:00 252.3 F 17.9  
28:00 252.0 F 17.7  
29:00 251.5 F 17.5  
30:00 252.4 F 17.9

Min 251.4 F 17.5  
Max 252.8 F 18.0

VENTING CHAMBER

DRYING  
Drying Start  
mm:ss Degrees PSI  
0:00 206.4 F 0.0  
5:00 198.7 F 0.0  
10:00 203.7 F 0.0  
15:00 207.7 F 0.0  
20:00 206.9 F 0.0  
25:00 199.6 F 0.0  
30:00 213.2 F 0.0  
Drying Complete

FILLING: 0:58  
HEATING: 22:19  
STERILIZING: 30:00  
VENTING: 3:06  
DRYING: 30:00  
TOTAL CYCLE: 01:26:24  
CYCLE COMPLETE

07/20/2021 11:58 AM  
ACCEPTED: X  
REJECTED:  
Operator: U01