
Midmark M11 - v1.0.0
Total Cycles: 629

Name: MIDMARK
Sterilizer ID: V2178198
01/26/2021 10:23 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 127.5 F 0.0
2:00 129.2 F 0.7
4:00 149.7 F 2.4
6:00 175.4 F 6.4
8:00 199.8 F 1.8
10:00 215.9 F 5.9
12:00 217.0 F 3.3
14:00 226.6 F 5.6
16:00 233.8 F 7.7
18:00 245.9 F 13.1

STERILIZING
mm:ss Degrees PSI
0:00 251.1 F 16.0
1:00 252.6 F 16.5
2:00 251.0 F 15.9
3:00 252.8 F 16.5
4:00 252.8 F 16.5
5:00 252.8 F 16.6
6:00 253.0 F 16.6
7:00 252.9 F 16.8
8:00 252.5 F 16.8
9:00 252.1 F 16.5
10:00 252.7 F 16.5
11:00 252.9 F 16.8
12:00 252.3 F 16.8
13:00 252.1 F 16.5
14:00 252.7 F 16.7
15:00 252.1 F 16.7
16:00 252.7 F 16.6
17:00 252.8 F 16.8
18:00 251.9 F 16.4
19:00 252.6 F 16.7
20:00 252.2 F 16.8
21:00 252.6 F 16.6
22:00 252.7 F 16.9
23:00 252.3 F 16.4
24:00 253.0 F 16.8
25:00 252.1 F 16.5
26:00 252.6 F 16.7
27:00 252.0 F 16.7
28:00 252.8 F 16.6
29:00 252.3 F 16.9
30:00 252.6 F 16.6

Min 250.9 F 15.9
Max 253.1 F 16.9

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 206.0 F 0.0
5:00 212.8 F 0.0
10:00 216.0 F 0.0
15:00 217.8 F 0.0
20:00 188.5 F 0.0
25:00 189.6 F 0.0
30:00 186.4 F 0.0
Drying Complete

FILLING: 0:58

HEATING: 19:36
STERILIZING: 30:00
VENTING: 2:53
DRYING: 30:00
TOTAL CYCLE: 01:23:26
CYCLE COMPLETE

01/26/2021 11:47 AM
ACCEPTED: X
REJECTED:
Operator: U01