
Midmark M11 - v1.0.0
Total Cycles: 602

Name: MIDMARK
Sterilizer ID: V2178198
01/13/2021 02:09 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING

mm:ss	Degrees	PSI
0:00	90.1 F	0.0
2:00	93.5 F	0.0
4:00	113.9 F	0.1
6:00	156.2 F	0.2
8:00	182.4 F	2.7
10:00	208.0 F	6.6
12:00	218.2 F	3.4
14:00	225.6 F	7.6
16:00	228.1 F	5.6
18:00	225.1 F	4.4
20:00	240.0 F	10.3
22:00	251.9 F	16.1

STERILIZING

mm:ss	Degrees	PSI
0:00	251.0 F	15.9
1:00	251.8 F	15.9
2:00	252.3 F	16.3
3:00	251.9 F	16.4
4:00	251.8 F	16.3
5:00	251.9 F	16.1
6:00	252.2 F	16.3
7:00	252.3 F	16.5
8:00	252.8 F	16.7
9:00	251.8 F	16.2
10:00	252.3 F	16.3
11:00	252.9 F	16.6
12:00	251.8 F	16.2
13:00	252.4 F	16.4
14:00	252.3 F	16.6
15:00	252.1 F	16.3
16:00	252.8 F	16.7
17:00	252.1 F	16.2
18:00	252.7 F	16.6
19:00	251.9 F	16.2
20:00	252.9 F	16.6
21:00	251.8 F	16.1
22:00	252.6 F	16.6
23:00	251.8 F	16.1
24:00	252.9 F	16.6
25:00	251.8 F	16.1
26:00	252.6 F	16.6
27:00	251.8 F	16.1
28:00	252.6 F	16.6
29:00	251.9 F	16.2
30:00	252.9 F	16.6

Min 251.0 F 15.8
Max 253.0 F 16.7

VENTING CHAMBER

DRYING

mm:ss	Degrees	PSI
0:00	174.3 F	0.0
5:00	233.0 F	0.0
10:00	217.6 F	0.0
15:00	200.3 F	0.0
20:00	168.1 F	0.0

Drying Start
Dry Cycle Aborted

FILLING: 0:56

HEATING: 22:28
STERILIZING: 30:00
VENTING: 2:53
DRYING: 22:60
TOTAL CYCLE: 01:19:17

01/13/2021 03:29 PM
ACCEPTED: X
REJECTED:
Operator: U01