
Midmark M11 - v1.0.0
Total Cycles: 582

Name: MIDMARK
Sterilizer ID: V2178198
12/22/2020 12:37 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 106.8 F 0.0
2:00 116.3 F 0.1
4:00 132.0 F 0.2
6:00 159.5 F 1.9
8:00 186.3 F 5.5
10:00 209.4 F 1.4
12:00 219.4 F 5.3
14:00 222.9 F 4.0
16:00 232.1 F 7.0
18:00 231.8 F 6.8
20:00 244.1 F 12.2

STERILIZING
mm:ss Degrees PSI
0:00 251.2 F 15.9
1:00 252.3 F 16.3
2:00 251.3 F 16.0
3:00 252.3 F 16.3
4:00 251.3 F 15.8
5:00 252.5 F 16.5
6:00 252.0 F 16.1
7:00 251.3 F 16.0
8:00 252.6 F 16.5
9:00 252.2 F 16.3
10:00 251.6 F 15.9
11:00 251.5 F 16.1
12:00 252.6 F 16.5
13:00 252.4 F 16.4
14:00 251.9 F 16.1
15:00 251.4 F 15.8
16:00 251.6 F 16.2
17:00 252.6 F 16.5
18:00 252.5 F 16.4
19:00 252.1 F 16.2
20:00 251.7 F 16.0
21:00 251.3 F 15.9
22:00 251.8 F 16.3
23:00 252.6 F 16.6
24:00 252.5 F 16.5
25:00 252.3 F 16.3
26:00 251.8 F 16.1
27:00 251.5 F 15.9
28:00 251.4 F 16.0
29:00 252.1 F 16.5
30:00 252.6 F 16.6

Min 251.2 F 15.8
Max 252.7 F 16.6

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 196.7 F 0.0
5:00 185.8 F 0.0
10:00 177.8 F 0.0
15:00 178.1 F 0.0
20:00 178.6 F 0.0
25:00 178.7 F 0.0
Dry Cycle Aborted

FILLING: 0:53

HEATING: 21:58
STERILIZING: 30:00
VENTING: 2:54
DRYING: 25:50
TOTAL CYCLE: 01:21:35

12/22/2020 01:58 PM
ACCEPTED: X
REJECTED:
Operator: U01