
Midmark M11 - v1.0.0
Total Cycles: 596

Name: MIDMARK
Sterilizer ID: V2178198
01/12/2021 08:28 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING

mm:ss	Degrees	PSI
0:00	72.9 F	0.0
2:00	77.7 F	0.0
4:00	117.0 F	0.1
6:00	133.1 F	0.2
8:00	166.6 F	0.3
10:00	190.3 F	2.2
12:00	205.1 F	5.6
14:00	212.3 F	2.6
16:00	224.0 F	5.5
18:00	223.9 F	3.9
20:00	231.1 F	6.3
22:00	234.1 F	7.6
24:00	247.2 F	13.5

STERILIZING

mm:ss	Degrees	PSI
0:00	251.5 F	15.6
1:00	252.0 F	15.7
2:00	252.7 F	16.1
3:00	252.1 F	16.1
4:00	251.8 F	15.6
5:00	252.8 F	16.2
6:00	251.6 F	15.6
7:00	252.9 F	16.2
8:00	251.7 F	15.6
9:00	252.8 F	16.3
10:00	252.2 F	15.8
11:00	252.0 F	16.0
12:00	252.8 F	16.1
13:00	251.8 F	15.6
14:00	252.6 F	16.3
15:00	252.7 F	16.1
16:00	251.8 F	15.6
17:00	252.6 F	16.3
18:00	252.7 F	16.1
19:00	251.9 F	15.7
20:00	252.2 F	16.1
21:00	252.9 F	16.2
22:00	252.1 F	15.8
23:00	251.7 F	15.8
24:00	253.0 F	16.3
25:00	252.4 F	16.0
26:00	251.7 F	15.6
27:00	252.6 F	16.3
28:00	252.8 F	16.2
29:00	252.1 F	15.8
30:00	251.6 F	15.7

Min 251.5 F 15.5
Max 253.0 F 16.3

VENTING CHAMBER

DRYING

Drying Start	mm:ss	Degrees	PSI
	0:00	174.2 F	0.0
	5:00	184.4 F	0.0
	10:00	185.9 F	0.0
	15:00	186.6 F	0.0
	20:00	188.5 F	0.0
	25:00	190.0 F	0.0
	30:00	170.1 F	0.0

Drying Complete

FILLING: 0:48
HEATING: 25:16
STERILIZING: 30:00
VENTING: 2:46
DRYING: 30:00
TOTAL CYCLE: 01:28:51
CYCLE COMPLETE

01/12/2021 09:56 AM
ACCEPTED: X
REJECTED:
Operator: U01