
Midmark M11 - v1.0.0
Total Cycles: 599

Name: MIDMARK
Sterilizer ID: V2178198
01/13/2021 07:18 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 72.8 F 0.0
2:00 75.8 F 0.0
4:00 99.6 F 0.0
6:00 131.5 F 0.1
8:00 151.8 F 0.3
10:00 179.3 F 2.4
12:00 198.8 F 5.2
14:00 210.2 F 2.0
16:00 220.0 F 5.0
18:00 224.5 F 4.7
20:00 230.1 F 6.5
22:00 227.7 F 5.2
24:00 241.0 F 10.6
26:00 252.3 F 16.2

STERILIZING
mm:ss Degrees PSI
0:00 251.4 F 15.8
1:00 252.2 F 16.1
2:00 252.3 F 16.4
3:00 251.6 F 15.7
4:00 252.6 F 16.3
5:00 251.5 F 15.7
6:00 252.6 F 16.3
7:00 251.7 F 15.8
8:00 252.6 F 16.4
9:00 252.1 F 16.0
10:00 251.7 F 16.0
11:00 252.6 F 16.3
12:00 251.8 F 15.9
13:00 252.3 F 16.4
14:00 252.4 F 16.2
15:00 251.6 F 15.8
16:00 252.5 F 16.4
17:00 252.4 F 16.2
18:00 251.7 F 15.8
19:00 252.3 F 16.4
20:00 252.5 F 16.3
21:00 251.8 F 15.9
22:00 251.8 F 16.1
23:00 252.6 F 16.4
24:00 252.1 F 16.1
25:00 251.5 F 15.8
26:00 252.7 F 16.5
27:00 252.3 F 16.2
28:00 251.8 F 15.9
29:00 251.8 F 16.1
30:00 252.7 F 16.4

Min 251.4 F 15.7
Max 252.8 F 16.5

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 193.3 F 0.0
5:00 193.4 F 0.0
10:00 177.1 F 0.0
15:00 174.9 F 0.0
20:00 172.7 F 0.0
Dry Cycle Aborted

FILLING: 0:53
HEATING: 26:27
STERILIZING: 30:00
VENTING: 2:50
DRYING: 21:33
TOTAL CYCLE: 01:21:42

01/13/2021 08:41 AM
ACCEPTED: X
REJECTED:
Operator: U01