
Midmark M11 - v1.0.0
Total Cycles: 593

Name: MIDMARK
Sterilizer ID: V2178198
01/11/2021 10:00 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 103.6 F 0.0
2:00 106.3 F 0.1
4:00 129.8 F 1.3
6:00 152.6 F 3.8
8:00 183.7 F 5.2
10:00 205.4 F 2.9
12:00 220.4 F 7.2
14:00 225.8 F 5.1
16:00 225.9 F 4.4
18:00 237.8 F 9.2
20:00 250.2 F 15.1

STERILIZING
mm:ss Degrees PSI
0:00 251.8 F 15.8
1:00 252.4 F 16.3
2:00 251.5 F 16.0
3:00 252.0 F 16.0
4:00 252.2 F 16.4
5:00 251.9 F 15.9
6:00 252.3 F 16.4
7:00 252.0 F 16.0
8:00 251.7 F 16.1
9:00 252.3 F 16.2
10:00 251.6 F 15.8
11:00 252.1 F 16.4
12:00 252.4 F 16.2
13:00 251.8 F 15.9
14:00 251.8 F 16.2
15:00 252.5 F 16.4
16:00 252.2 F 16.1
17:00 251.5 F 15.8
18:00 251.8 F 16.2
19:00 252.6 F 16.4
20:00 252.2 F 16.2
21:00 251.7 F 15.9
22:00 251.4 F 16.0
23:00 252.4 F 16.4
24:00 252.3 F 16.3
25:00 252.0 F 16.0
26:00 251.5 F 15.8
27:00 251.7 F 16.1
28:00 252.5 F 16.4
29:00 252.3 F 16.3
30:00 252.1 F 16.0

Min 251.2 F 15.7
Max 252.6 F 16.5

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 184.9 F 0.0
5:00 181.3 F 0.0
10:00 175.7 F 0.0
15:00 176.6 F 0.0
20:00 178.3 F 0.0
25:00 183.5 F 0.0
30:00 181.1 F 0.0

Drying Complete

FILLING: 0:51
HEATING: 20:41
STERILIZING: 30:00
VENTING: 2:47
DRYING: 30:00
TOTAL CYCLE: 01:24:18
CYCLE COMPLETE

01/11/2021 11:32 AM
ACCEPTED: X
REJECTED:
Operator: U01