
Midmark M11 - v1.0.5
Total Cycles: 15

Name: RITTER2
Sterilizer ID: V2313228
04/26/2021 08:20 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 118.9 F 0.0
2:00 121.5 F 0.8
4:00 129.7 F 3.2
6:00 155.2 F 6.6
8:00 190.1 F 1.1
10:00 206.6 F 5.2
12:00 218.3 F 3.8
14:00 228.6 F 7.3
16:00 230.9 F 6.9
18:00 243.3 F 12.3
20:00 252.0 F 16.5

STERILIZING
mm:ss Degrees PSI
0:00 251.4 F 16.5
1:00 252.9 F 17.1
2:00 252.0 F 16.6
3:00 251.6 F 16.6
4:00 252.7 F 17.1
5:00 252.5 F 16.9
6:00 252.1 F 16.6
7:00 251.7 F 16.4
8:00 251.5 F 16.5
9:00 252.0 F 16.9
10:00 252.6 F 17.2
11:00 252.8 F 17.1
12:00 252.8 F 17.1
13:00 252.7 F 17.0
14:00 252.7 F 17.0
15:00 252.5 F 16.9
16:00 252.5 F 16.9
17:00 252.4 F 16.8
18:00 252.4 F 16.8
19:00 252.4 F 16.8
20:00 252.4 F 16.8
21:00 252.4 F 16.8
22:00 252.3 F 16.8
23:00 252.3 F 16.8
24:00 252.3 F 16.8
25:00 252.3 F 16.8
26:00 252.3 F 16.8
27:00 252.4 F 16.8
28:00 252.4 F 16.8
29:00 252.4 F 16.8
30:00 252.4 F 16.9

Min 251.4 F 16.3
Max 253.2 F 17.2

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 192.8 F 0.0
5:00 193.8 F 0.0
10:00 190.7 F 0.0
15:00 191.8 F 0.0
20:00 189.9 F 0.0
25:00 162.3 F 0.0
30:00 163.1 F 0.0

Drying Complete

FILLING: 1:00
HEATING: 20:10
STERILIZING: 30:00
VENTING: 2:58
DRYING: 30:00
TOTAL CYCLE: 01:24:08
CYCLE COMPLETE

04/26/2021 09:53 AM
ACCEPTED: X
REJECTED:
Operator: U01