
Midmark M11 - v1.0.5
Total Cycles: 59

Name: RITTER2
Sterilizer ID: V2313228
05/31/2021 01:46 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 164.3 F 0.0
2:00 168.1 F 2.9
4:00 196.6 F 4.1
6:00 218.1 F 4.9
8:00 223.9 F 4.4
10:00 228.4 F 5.5
12:00 239.1 F 10.0
14:00 252.0 F 16.3

STERILIZING
mm:ss Degrees PSI
0:00 251.6 F 16.0
1:00 252.0 F 16.1
2:00 252.6 F 16.4
3:00 252.0 F 16.5
4:00 252.1 F 16.2
5:00 252.3 F 16.6
6:00 252.4 F 16.3
7:00 251.6 F 16.0
8:00 252.8 F 16.6
9:00 252.3 F 16.3
10:00 251.7 F 16.0
11:00 252.3 F 16.7
12:00 252.7 F 16.6
13:00 252.3 F 16.3
14:00 251.9 F 16.1
15:00 251.6 F 16.1
16:00 252.3 F 16.7
17:00 252.8 F 16.6
18:00 252.4 F 16.5
19:00 252.2 F 16.3
20:00 251.9 F 16.1
21:00 251.6 F 16.0
22:00 251.7 F 16.3
23:00 252.4 F 16.7
24:00 252.8 F 16.6
25:00 252.6 F 16.5
26:00 252.4 F 16.4
27:00 252.1 F 16.2
28:00 252.0 F 16.1
29:00 251.6 F 16.0
30:00 251.6 F 16.1

Min 251.4 F 15.9
Max 252.9 F 16.7

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 180.9 F 0.0
5:00 204.0 F 0.0
10:00 176.2 F 0.0
15:00 173.1 F 0.0
Dry Cycle Aborted

FILLING: 2:24
HEATING: 14:28
STERILIZING: 30:00
VENTING: 2:50
DRYING: 17:47
TOTAL CYCLE: 01:07:29

05/31/2021 02:53 PM

ACCEPTED: X

REJECTED:

Operator: U01