
Midmark M11 - v1.0.0
Total Cycles: 594

Name: MIDMARK
Sterilizer ID: V2178198
01/11/2021 11:54 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 137.1 F 0.0
2:00 151.2 F 1.1
4:00 167.5 F 4.2
6:00 197.3 F 1.8
8:00 213.1 F 4.8
10:00 220.6 F 3.6
12:00 230.0 F 6.3
14:00 233.8 F 7.6
16:00 247.1 F 13.6

STERILIZING
mm:ss Degrees PSI
0:00 251.4 F 15.9
1:00 252.3 F 16.2
2:00 252.1 F 16.4
3:00 252.0 F 16.0
4:00 252.3 F 16.5
5:00 252.2 F 16.1
6:00 251.8 F 16.3
7:00 252.5 F 16.3
8:00 251.6 F 15.9
9:00 252.3 F 16.5
10:00 252.6 F 16.3
11:00 251.9 F 16.0
12:00 251.4 F 16.0
13:00 252.7 F 16.5
14:00 252.6 F 16.3
15:00 252.2 F 16.1
16:00 251.5 F 15.8
17:00 251.9 F 16.3
18:00 252.7 F 16.5
19:00 252.5 F 16.3
20:00 252.1 F 16.1
21:00 251.5 F 15.9
22:00 251.7 F 16.1
23:00 252.8 F 16.6
24:00 252.6 F 16.4
25:00 252.4 F 16.2
26:00 252.0 F 16.0
27:00 251.4 F 15.8
28:00 251.8 F 16.2
29:00 252.7 F 16.5
30:00 252.7 F 16.4

Min 251.3 F 15.8
Max 252.8 F 16.6

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 202.4 F 0.0
5:00 195.2 F 0.0
10:00 186.0 F 0.0
15:00 186.0 F 0.0
20:00 186.3 F 0.0
25:00 191.6 F 0.0
30:00 188.7 F 0.0
Drying Complete

FILLING: 0:56
HEATING: 17:20

STERILIZING: 30:00
VENTING: 2:49
DRYING: 30:00
TOTAL CYCLE: 01:21:05
CYCLE COMPLETE

01/11/2021 01:15 PM
ACCEPTED: X
REJECTED:
Operator: U01