
Midmark M11 - v1.0.5
Total Cycles: 98

Name: RITTER2
Sterilizer ID: V2313228
06/14/2021 03:37 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 153.3 F 0.0
2:00 155.8 F 2.8
4:00 176.7 F 6.8
6:00 205.2 F 3.3
8:00 223.7 F 5.8
10:00 228.1 F 7.0
12:00 229.9 F 6.6
14:00 243.5 F 12.4

STERILIZING
mm:ss Degrees PSI
0:00 251.5 F 16.6
1:00 252.8 F 17.3
2:00 252.7 F 17.2
3:00 252.0 F 17.1
4:00 252.0 F 17.2
5:00 252.5 F 17.3
6:00 252.9 F 17.4
7:00 252.6 F 17.6
8:00 251.6 F 17.2
9:00 252.6 F 17.3
10:00 252.4 F 17.6
11:00 251.8 F 17.2
12:00 252.6 F 17.5
13:00 251.7 F 17.1
14:00 252.5 F 17.5
15:00 251.8 F 17.1
16:00 252.5 F 17.6
17:00 252.1 F 17.3
18:00 252.1 F 17.6
19:00 252.5 F 17.5
20:00 251.7 F 17.1
21:00 252.6 F 17.6
22:00 252.6 F 17.4
23:00 251.7 F 17.1
24:00 252.5 F 17.6
25:00 252.5 F 17.4
26:00 251.7 F 17.1
27:00 252.6 F 17.6
28:00 252.5 F 17.4
29:00 251.8 F 17.2
30:00 252.3 F 17.6

Min 250.8 F 16.5
Max 253.4 F 17.7

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 180.5 F 0.0
5:00 202.4 F 0.0
Dry Cycle Aborted

FILLING: 1:45
HEATING: 15:50
STERILIZING: 30:00
VENTING: 2:57
DRYING: 8:40
TOTAL CYCLE: 00:59:12

06/14/2021 04:36 PM

ACCEPTED: X
REJECTED:
Operator: U01