
Midmark M11 - v1.0.0
Total Cycles: 672

Name: MIDMARK
Sterilizer ID: V2178198
02/17/2021 12:09 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 179.2 F 0.0
2:00 174.9 F 3.5
4:00 186.1 F 7.4
6:00 208.1 F 4.1
8:00 225.4 F 4.6
10:00 229.0 F 7.2
12:00 231.4 F 6.3
14:00 244.6 F 12.2
16:00 252.2 F 16.2

STERILIZING
mm:ss Degrees PSI
0:00 251.3 F 15.9
1:00 251.4 F 15.9
2:00 251.7 F 15.9
3:00 252.5 F 16.1
4:00 252.6 F 16.3
5:00 251.1 F 16.0
6:00 252.0 F 16.1
7:00 252.3 F 16.4
8:00 251.6 F 16.0
9:00 253.1 F 16.7
10:00 251.9 F 16.3
11:00 252.8 F 16.5
12:00 253.1 F 16.7
13:00 252.0 F 16.3
14:00 252.7 F 16.5
15:00 252.8 F 16.7
16:00 251.9 F 16.3
17:00 253.3 F 16.7
18:00 251.9 F 16.4
19:00 252.3 F 16.5
20:00 252.6 F 16.7
21:00 252.0 F 16.4
22:00 253.2 F 16.7
23:00 251.9 F 16.3
24:00 252.4 F 16.7
25:00 251.8 F 16.3
26:00 253.0 F 16.6
27:00 251.8 F 16.5
28:00 252.6 F 16.5
29:00 252.2 F 16.7
30:00 252.6 F 16.5

Min 250.8 F 15.9
Max 253.5 F 16.7

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 197.1 F 0.0
5:00 235.2 F 0.0
10:00 224.3 F 0.0
15:00 232.1 F 0.0
20:00 217.3 F 0.0
Dry Cycle Aborted

FILLING: 1:14
HEATING: 16:25
STERILIZING: 30:00
VENTING: 2:53

DRYING: 20:39
TOTAL CYCLE: 01:11:11

02/17/2021 01:20 PM
ACCEPTED: X
REJECTED:
Operator: U01