

-----  
Midmark M11 - v1.0.0  
Total Cycles: 598

Name: MIDMARK  
Sterilizer ID: V2178198  
01/12/2021 01:22 PM

BEGIN PACKS CYCLE  
Temp: 250 Degrees F  
Time: 30 Minutes  
Dry: 30 Minutes

FILLING CHAMBER

HEATING  
mm:ss Degrees PSI  
0:00 97.9 F 0.0  
2:00 97.8 F 0.0  
4:00 109.1 F 0.2  
6:00 147.1 F 1.4  
8:00 174.6 F 4.9  
10:00 212.8 F 2.5  
12:00 215.1 F 5.2  
14:00 222.6 F 3.7  
16:00 231.1 F 7.4  
18:00 230.6 F 6.5  
20:00 243.8 F 12.1

STERILIZING  
mm:ss Degrees PSI  
0:00 251.4 F 16.0  
1:00 252.2 F 16.2  
2:00 252.6 F 16.5  
3:00 251.5 F 15.9  
4:00 252.5 F 16.4  
5:00 251.5 F 16.0  
6:00 252.6 F 16.4  
7:00 251.5 F 15.9  
8:00 252.8 F 16.5  
9:00 252.0 F 16.1  
10:00 252.0 F 16.4  
11:00 252.6 F 16.4  
12:00 251.9 F 16.0  
13:00 251.9 F 16.3  
14:00 252.6 F 16.5  
15:00 252.2 F 16.2  
16:00 251.5 F 15.9  
17:00 252.6 F 16.6  
18:00 252.7 F 16.4  
19:00 252.0 F 16.1  
20:00 251.5 F 16.0  
21:00 252.6 F 16.6  
22:00 252.7 F 16.5  
23:00 252.2 F 16.2  
24:00 251.7 F 16.0  
25:00 251.9 F 16.3  
26:00 252.8 F 16.6  
27:00 252.2 F 16.3  
28:00 252.0 F 16.1  
29:00 251.5 F 15.9  
30:00 252.5 F 16.6

Min 251.4 F 15.9  
Max 252.9 F 16.6

VENTING CHAMBER

DRYING  
Drying Start  
mm:ss Degrees PSI  
0:00 196.3 F 0.0  
5:00 208.6 F 0.0  
10:00 201.7 F 0.0  
15:00 182.5 F 0.0  
20:00 179.6 F 0.0  
25:00 178.6 F 0.0  
30:00 178.1 F 0.0  
Drying Complete

FILLING: 0:51  
HEATING: 21:55  
STERILIZING: 30:00  
VENTING: 2:53  
DRYING: 30:00  
TOTAL CYCLE: 01:25:38  
CYCLE COMPLETE

01/12/2021 02:57 PM  
ACCEPTED: X  
REJECTED:  
Operator: U01