

-----  
Midmark M11 - v1.0.0  
Total Cycles: 589

Name: MIDMARK  
Sterilizer ID: V2178198  
12/30/2020 09:30 AM

BEGIN PACKS CYCLE  
Temp: 250 Degrees F  
Time: 30 Minutes  
Dry: 30 Minutes

FILLING CHAMBER

HEATING  
mm:ss Degrees PSI  
0:00 146.0 F 0.0  
2:00 164.8 F 1.0  
4:00 193.2 F 3.6  
6:00 210.2 F 4.8  
8:00 214.4 F 3.6  
10:00 226.8 F 5.3  
12:00 230.5 F 6.5  
14:00 228.4 F 5.3  
16:00 242.0 F 10.9  
18:00 252.6 F 16.1

STERILIZING  
mm:ss Degrees PSI  
0:00 251.5 F 16.2  
1:00 252.5 F 16.4  
2:00 252.0 F 16.1  
3:00 251.6 F 16.4  
4:00 252.3 F 16.4  
5:00 252.5 F 16.4  
6:00 252.6 F 16.4  
7:00 252.3 F 16.4  
8:00 252.8 F 16.4  
9:00 252.5 F 16.4  
10:00 252.0 F 16.5  
11:00 251.1 F 16.3  
12:00 251.3 F 16.1  
13:00 253.3 F 17.0  
14:00 252.4 F 16.6  
15:00 251.7 F 16.6  
16:00 253.3 F 17.0  
17:00 252.9 F 16.8  
18:00 251.8 F 16.5  
19:00 252.1 F 16.9  
20:00 253.5 F 17.0  
21:00 252.9 F 16.8  
22:00 251.9 F 16.5  
23:00 252.0 F 16.9  
24:00 253.3 F 17.0  
25:00 252.7 F 16.8  
26:00 251.8 F 16.6  
27:00 251.6 F 16.7  
28:00 253.3 F 17.0  
29:00 253.0 F 16.9  
30:00 252.3 F 16.6

Min 251.0 F 16.1  
Max 253.8 F 17.1

VENTING CHAMBER

DRYING  
Drying Start  
mm:ss Degrees PSI  
0:00 198.9 F 0.0  
5:00 193.1 F 0.0  
10:00 192.5 F 0.0  
15:00 200.9 F 0.0  
20:00 204.2 F 0.0  
Dry Cycle Aborted

FILLING: 0:58  
HEATING: 18:48  
STERILIZING: 30:00

VENTING: 2:46  
DRYING: 23:08  
TOTAL CYCLE: 01:15:40

12/30/2020 10:45 AM  
ACCEPTED: X  
REJECTED:  
Operator: U01