
Midmark M11 - v1.0.5
Total Cycles: 60

Name: RITTER2
Sterilizer ID: V2313228
05/31/2021 02:54 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 175.0 F 0.0
2:00 175.4 F 3.6
4:00 203.5 F 3.2
6:00 222.9 F 6.3
8:00 227.0 F 5.7
10:00 228.8 F 5.8
12:00 244.6 F 12.6

STERILIZING
mm:ss Degrees PSI
0:00 251.5 F 16.1
1:00 252.2 F 16.3
2:00 252.4 F 16.7
3:00 252.0 F 16.2
4:00 251.7 F 16.5
5:00 252.3 F 16.5
6:00 251.8 F 16.1
7:00 251.5 F 16.3
8:00 252.5 F 16.7
9:00 252.1 F 16.5
10:00 251.8 F 16.2
11:00 251.5 F 16.0
12:00 251.4 F 16.3
13:00 252.2 F 16.8
14:00 252.5 F 16.7
15:00 252.4 F 16.6
16:00 252.2 F 16.5
17:00 251.9 F 16.3
18:00 251.6 F 16.2
19:00 251.5 F 16.1
20:00 251.3 F 16.0
21:00 251.3 F 16.2
22:00 251.5 F 16.5
23:00 252.1 F 16.8
24:00 252.3 F 16.8
25:00 252.5 F 16.7
26:00 252.4 F 16.7
27:00 252.3 F 16.6
28:00 252.4 F 16.6
29:00 252.3 F 16.6
30:00 252.3 F 16.5

Min 251.1 F 16.0
Max 252.8 F 16.8

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 202.7 F 0.0
5:00 190.7 F 0.0
Dry Cycle Aborted

FILLING: 2:40
HEATING: 13:39
STERILIZING: 30:00
VENTING: 2:49
DRYING: 5:08
TOTAL CYCLE: 00:54:16

05/31/2021 03:48 PM
ACCEPTED: X

REJECTED:
Operator: U01