
Midmark M11 - v1.0.5
Total Cycles: 6

Name: RITTER2
Sterilizer ID: V2313228
04/21/2021 02:42 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 166.9 F 0.0
2:00 170.7 F 3.1
4:00 200.3 F 3.9
6:00 219.3 F 5.8
8:00 224.6 F 5.2
10:00 227.8 F 5.6
12:00 242.9 F 12.2

STERILIZING
mm:ss Degrees PSI
0:00 251.7 F 16.6
1:00 252.3 F 16.7
2:00 252.8 F 17.1
3:00 251.7 F 16.5
4:00 252.7 F 17.0
5:00 251.6 F 16.5
6:00 252.9 F 17.1
7:00 252.0 F 16.6
8:00 252.5 F 17.2
9:00 252.3 F 16.9
10:00 251.6 F 16.5
11:00 252.8 F 17.2
12:00 252.4 F 16.9
13:00 251.7 F 16.5
14:00 252.8 F 17.2
15:00 252.6 F 17.0
16:00 252.0 F 16.6
17:00 251.9 F 16.8
18:00 252.8 F 17.1
19:00 252.4 F 16.9
20:00 251.9 F 16.6
21:00 252.1 F 17.0
22:00 252.8 F 17.1
23:00 252.4 F 16.9
24:00 251.9 F 16.6
25:00 252.0 F 17.0
26:00 252.8 F 17.1
27:00 252.4 F 16.9
28:00 251.9 F 16.6
29:00 251.7 F 16.8
30:00 252.7 F 17.2

Min 251.5 F 16.5
Max 253.2 F 17.2

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 178.9 F 0.0
5:00 192.3 F 0.0
10:00 191.7 F 0.0
15:00 193.4 F 0.0
20:00 195.1 F 0.0
25:00 199.1 F 0.0
30:00 171.6 F 0.0
Drying Complete

FILLING: 2:25
HEATING: 13:51
STERILIZING: 30:00
VENTING: 2:55

DRYING: 30:00
TOTAL CYCLE: 01:19:11
CYCLE COMPLETE

04/21/2021 04:01 PM
ACCEPTED: X
REJECTED:
Operator: U01