

-----  
Midmark M11 - v1.0.0  
Total Cycles: 673

Name: MIDMARK  
Sterilizer ID: V2178198  
02/17/2021 01:21 PM

BEGIN PACKS CYCLE  
Temp: 250 Degrees F  
Time: 30 Minutes  
Dry: 30 Minutes

FILLING CHAMBER

HEATING  
mm:ss Degrees PSI  
0:00 183.9 F 0.0  
2:00 179.4 F 3.4  
4:00 204.6 F 5.0  
6:00 223.7 F 5.2  
8:00 222.4 F 3.9  
10:00 227.3 F 5.2  
12:00 238.3 F 9.2  
14:00 251.1 F 15.3

STERILIZING  
mm:ss Degrees PSI  
0:00 252.3 F 15.9  
1:00 251.9 F 15.8  
2:00 251.5 F 15.8  
3:00 251.9 F 15.9  
4:00 252.6 F 16.0  
5:00 252.6 F 16.2  
6:00 252.0 F 16.2  
7:00 251.6 F 15.9  
8:00 252.5 F 16.2  
9:00 251.3 F 16.0  
10:00 252.1 F 16.1  
11:00 251.8 F 16.2  
12:00 253.1 F 16.5  
13:00 251.6 F 16.2  
14:00 253.0 F 16.6  
15:00 252.2 F 16.3  
16:00 252.5 F 16.6  
17:00 252.9 F 16.5  
18:00 251.9 F 16.2  
19:00 253.0 F 16.7  
20:00 253.1 F 16.5  
21:00 251.8 F 16.2  
22:00 253.2 F 16.7  
23:00 252.2 F 16.4  
24:00 251.9 F 16.2  
25:00 253.0 F 16.7  
26:00 252.5 F 16.4  
27:00 251.9 F 16.2  
28:00 252.6 F 16.7  
29:00 252.9 F 16.4  
30:00 251.8 F 16.2

Min 251.0 F 15.8  
Max 253.6 F 16.7

VENTING CHAMBER

DRYING  
Drying Start  
mm:ss Degrees PSI  
0:00 197.2 F 0.0  
5:00 240.0 F 0.0  
10:00 227.1 F 0.0  
15:00 235.0 F 0.0  
20:00 201.4 F 0.0  
Dry Cycle Aborted

FILLING: 2:00  
HEATING: 15:05  
STERILIZING: 30:00  
VENTING: 2:52  
DRYING: 21:01

TOTAL CYCLE: 01:10:59

02/17/2021 02:32 PM

ACCEPTED: X

REJECTED:

Operator: U01