
Midmark M11 - v1.0.5
Total Cycles: 5

Name: RITTER2
Sterilizer ID: V2313228
04/21/2021 01:27 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 163.9 F 0.0
2:00 168.4 F 2.9
4:00 195.1 F 4.6
6:00 216.7 F 5.1
8:00 224.0 F 4.8
10:00 227.9 F 5.6
12:00 241.9 F 11.7

STERILIZING
mm:ss Degrees PSI
0:00 252.0 F 16.5
1:00 252.0 F 16.6
2:00 252.4 F 16.8
3:00 252.8 F 17.0
4:00 251.4 F 16.6
5:00 252.0 F 16.6
6:00 252.6 F 17.0
7:00 251.5 F 16.5
8:00 252.6 F 16.9
9:00 251.5 F 16.4
10:00 252.4 F 17.0
11:00 251.7 F 16.5
12:00 252.6 F 17.0
13:00 251.9 F 16.6
14:00 252.1 F 17.1
15:00 252.1 F 16.7
16:00 251.4 F 16.7
17:00 252.4 F 16.9
18:00 251.7 F 16.5
19:00 252.6 F 17.1
20:00 252.1 F 16.7
21:00 251.3 F 16.6
22:00 252.4 F 17.0
23:00 251.8 F 16.6
24:00 251.5 F 16.8
25:00 252.4 F 16.9
26:00 251.8 F 16.5
27:00 252.3 F 17.1
28:00 252.3 F 16.8
29:00 251.4 F 16.4
30:00 252.4 F 17.1

Min 251.2 F 16.4
Max 252.9 F 17.1

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 179.3 F 0.0
5:00 190.1 F 0.0
10:00 158.6 F 0.0
15:00 155.9 F 0.0
20:00 155.7 F 0.0
25:00 146.3 F 0.0
Dry Cycle Aborted

FILLING: 2:25
HEATING: 13:50
STERILIZING: 30:00
VENTING: 2:54
DRYING: 25:21

TOTAL CYCLE: 01:14:30

04/21/2021 02:41 PM

ACCEPTED: X

REJECTED:

Operator: U01