

-----  
Midmark M11 - v1.0.5  
Total Cycles: 68

Name: RITTER2  
Sterilizer ID: V2313228  
06/07/2021 08:55 AM

BEGIN PACKS CYCLE  
Temp: 250 Degrees F  
Time: 30 Minutes  
Dry: 30 Minutes

FILLING CHAMBER

HEATING  
mm:ss Degrees PSI  
0:00 178.2 F 0.0  
2:00 176.6 F 1.6  
4:00 179.7 F 5.7  
6:00 204.7 F 2.7  
8:00 220.7 F 6.1  
10:00 227.6 F 7.0  
12:00 230.1 F 6.7  
14:00 243.6 F 12.7

STERILIZING  
mm:ss Degrees PSI  
0:00 251.2 F 16.6  
1:00 252.7 F 17.2  
2:00 251.3 F 16.5  
3:00 252.5 F 17.4  
4:00 252.1 F 16.9  
5:00 251.3 F 16.5  
6:00 252.0 F 17.2  
7:00 253.0 F 17.2  
8:00 252.4 F 17.0  
9:00 251.9 F 16.8  
10:00 251.4 F 16.6  
11:00 251.3 F 16.8  
12:00 252.0 F 17.2  
13:00 252.4 F 17.4  
14:00 253.0 F 17.3  
15:00 252.7 F 17.3  
16:00 252.9 F 17.3  
17:00 252.7 F 17.2  
18:00 253.0 F 17.2  
19:00 253.0 F 17.2  
20:00 252.7 F 17.2  
21:00 252.8 F 17.2  
22:00 252.7 F 17.2  
23:00 252.7 F 17.2  
24:00 252.9 F 17.3  
25:00 253.1 F 17.3  
26:00 252.8 F 17.4  
27:00 252.3 F 17.4  
28:00 252.1 F 17.2  
29:00 251.6 F 17.0  
30:00 251.3 F 16.7

Min 251.1 F 16.5  
Max 253.4 F 17.4

VENTING CHAMBER

DRYING  
Drying Start  
mm:ss Degrees PSI  
0:00 203.7 F 0.0  
5:00 207.0 F 0.0  
10:00 207.3 F 0.0  
15:00 205.5 F 0.0  
Dry Cycle Aborted

FILLING: 1:59  
HEATING: 15:53  
STERILIZING: 30:00  
VENTING: 2:52  
DRYING: 18:41  
TOTAL CYCLE: 01:09:25

06/07/2021 10:04 AM  
ACCEPTED:  
REJECTED:  
Reason: \_\_\_\_\_  
Operator: \_\_\_\_\_