
Midmark M11 - v1.0.0
Total Cycles: 580

Name: MIDMARK
Sterilizer ID: V2178198
12/22/2020 07:27 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 73.8 F 0.0
2:00 75.2 F 0.0
4:00 91.0 F 0.1
6:00 123.6 F 0.1
8:00 159.5 F 0.3
10:00 185.2 F 1.7
12:00 199.6 F 5.0
14:00 211.3 F 3.7
16:00 221.0 F 4.3
18:00 228.1 F 6.1
20:00 230.0 F 6.3
22:00 227.5 F 5.2
24:00 241.5 F 10.9
26:00 252.7 F 16.3

STERILIZING
mm:ss Degrees PSI
0:00 251.5 F 15.8
1:00 252.1 F 16.0
2:00 252.8 F 16.4
3:00 251.6 F 15.9
4:00 252.4 F 16.1
5:00 252.2 F 16.2
6:00 252.2 F 16.0
7:00 252.2 F 16.2
8:00 252.5 F 16.2
9:00 251.6 F 15.8
10:00 252.8 F 16.4
11:00 252.0 F 16.0
12:00 251.9 F 16.1
13:00 252.8 F 16.4
14:00 252.1 F 16.0
15:00 251.6 F 15.9
16:00 252.9 F 16.5
17:00 252.5 F 16.2
18:00 251.9 F 15.9
19:00 252.0 F 16.2
20:00 252.9 F 16.4
21:00 252.4 F 16.1
22:00 251.8 F 15.8
23:00 252.1 F 16.2
24:00 252.9 F 16.5
25:00 252.5 F 16.2
26:00 252.0 F 16.0
27:00 251.5 F 15.8
28:00 252.7 F 16.5
29:00 252.8 F 16.4
30:00 252.3 F 16.1

Min 251.5 F 15.7
Max 252.9 F 16.5

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 176.3 F 0.0
5:00 185.3 F 0.0
10:00 182.5 F 0.0
15:00 178.2 F 0.0
20:00 163.9 F 0.0
25:00 158.9 F 0.0

30:00 156.2 F 0.0
Drying Complete

FILLING: 0:49
HEATING: 26:17
STERILIZING: 30:00
VENTING: 2:48
DRYING: 30:00
TOTAL CYCLE: 01:29:55
CYCLE COMPLETE

12/22/2020 08:56 AM
ACCEPTED: X
REJECTED:
Operator: U01