
Midmark M11 - v1.0.0
Total Cycles: 577

Name: MIDMARK
Sterilizer ID: V2178198
12/21/2020 11:01 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 174.1 F 0.0
2:00 168.7 F 2.3
4:00 175.6 F 6.2
6:00 199.9 F 1.6
8:00 217.8 F 6.6
10:00 225.4 F 4.9
12:00 226.6 F 4.6
14:00 239.4 F 10.0
16:00 251.8 F 16.0

STERILIZING
mm:ss Degrees PSI
0:00 251.1 F 15.9
1:00 252.6 F 16.3
2:00 251.6 F 16.1
3:00 252.7 F 16.4
4:00 252.1 F 16.5
5:00 252.8 F 16.4
6:00 251.8 F 16.4
7:00 252.7 F 16.5
8:00 251.8 F 16.1
9:00 252.6 F 16.7
10:00 252.7 F 16.4
11:00 251.5 F 16.1
12:00 252.9 F 16.7
13:00 252.7 F 16.5
14:00 251.8 F 16.1
15:00 252.1 F 16.6
16:00 252.8 F 16.6
17:00 252.6 F 16.4
18:00 251.6 F 16.1
19:00 252.2 F 16.6
20:00 253.0 F 16.6
21:00 252.5 F 16.4
22:00 251.6 F 16.1
23:00 252.0 F 16.5
24:00 252.9 F 16.7
25:00 252.6 F 16.5
26:00 252.0 F 16.2
27:00 251.8 F 16.3
28:00 253.0 F 16.7
29:00 252.6 F 16.5
30:00 252.3 F 16.3

Min 251.0 F 15.8
Max 253.1 F 16.7

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 192.4 F 0.0
5:00 214.6 F 0.0
10:00 200.1 F 0.0
15:00 190.7 F 0.0
20:00 192.9 F 0.0
25:00 221.9 F 0.0
30:00 231.0 F 0.0
Drying Complete

FILLING: 0:59
HEATING: 16:34

STERILIZING: 30:00
VENTING: 2:53
DRYING: 30:00
TOTAL CYCLE: 01:20:26
CYCLE COMPLETE

12/21/2020 12:21 PM
ACCEPTED:
REJECTED:
Reason: _____
Operator: _____