
Midmark M11 - v1.0.5
Total Cycles: 55

Name: RITTER2
Sterilizer ID: V2313228
05/31/2021 08:16 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 119.4 F 0.0
2:00 133.3 F 1.0
4:00 167.1 F 3.3
6:00 191.3 F 6.8
8:00 185.9 F 1.5
10:00 204.4 F 5.9
12:00 214.7 F 2.9
14:00 226.5 F 6.4
16:00 231.9 F 7.1
18:00 244.2 F 12.5

STERILIZING
mm:ss Degrees PSI
0:00 251.8 F 16.5
1:00 252.1 F 16.4
2:00 252.6 F 16.6
3:00 252.7 F 16.9
4:00 252.0 F 16.6
5:00 252.4 F 16.6
6:00 252.3 F 16.9
7:00 252.1 F 16.7
8:00 252.0 F 16.7
9:00 252.2 F 16.8
10:00 251.5 F 16.4
11:00 252.4 F 17.0
12:00 252.1 F 16.8
13:00 251.8 F 16.4
14:00 251.9 F 16.9
15:00 252.0 F 16.9
16:00 252.4 F 17.1
17:00 253.0 F 16.9
18:00 252.5 F 16.7
19:00 252.3 F 16.9
20:00 252.9 F 17.4
21:00 253.4 F 17.3
22:00 252.5 F 17.3
23:00 252.4 F 17.2
24:00 253.2 F 17.0
25:00 252.9 F 16.9
26:00 252.4 F 16.8
27:00 252.2 F 16.7
28:00 251.5 F 16.9
29:00 252.9 F 17.2
30:00 253.3 F 17.4

Min 250.9 F 16.4
Max 253.9 F 17.4

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 189.6 F 0.0
5:00 186.0 F 0.0
10:00 202.6 F 0.0
15:00 203.6 F 0.0
20:00 210.8 F 0.0
25:00 218.4 F 0.0
30:00 203.8 F 0.0
Drying Complete

FILLING: 0:59

HEATING: 19:54
STERILIZING: 30:00
VENTING: 2:51
DRYING: 30:00
TOTAL CYCLE: 01:23:44
CYCLE COMPLETE

05/31/2021 09:41 AM
ACCEPTED: X
REJECTED:
Operator: U01