
Midmark M11 - v1.0.5
Total Cycles: 67

Name: RITTER2
Sterilizer ID: V2313228
06/07/2021 07:29 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 70.6 F 0.0
2:00 72.2 F 0.2
4:00 100.5 F 0.9
6:00 118.3 F 3.0
8:00 135.5 F 5.8
10:00 189.8 F 0.8
12:00 192.1 F 3.4
14:00 210.6 F 6.8
16:00 216.8 F 3.5
18:00 223.1 F 7.4
20:00 226.1 F 5.7
22:00 237.4 F 10.5
24:00 248.3 F 15.7

STERILIZING
mm:ss Degrees PSI
0:00 251.4 F 17.5
1:00 252.2 F 17.6
2:00 252.5 F 18.1
3:00 251.5 F 17.4
4:00 252.4 F 18.0
5:00 251.5 F 17.4
6:00 252.3 F 18.1
7:00 252.1 F 17.7
8:00 251.4 F 17.6
9:00 252.5 F 18.0
10:00 252.4 F 17.7
11:00 251.5 F 17.4
12:00 251.5 F 17.8
13:00 252.9 F 18.0
14:00 252.4 F 17.9
15:00 252.2 F 17.7
16:00 251.6 F 17.5
17:00 251.4 F 17.5
18:00 251.6 F 17.8
19:00 252.4 F 18.1
20:00 252.7 F 18.0
21:00 252.7 F 17.9
22:00 252.6 F 17.8
23:00 252.2 F 17.8
24:00 252.0 F 17.7
25:00 251.8 F 17.6
26:00 251.8 F 17.5
27:00 251.6 F 17.5
28:00 251.5 F 17.4
29:00 251.4 F 17.5
30:00 251.5 F 17.6

Min 251.2 F 17.4
Max 253.4 F 18.1

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 205.7 F 0.0
5:00 191.9 F 0.0
10:00 199.4 F 0.0
15:00 202.8 F 0.0
20:00 171.5 F 0.0
25:00 160.5 F 0.0
Dry Cycle Aborted

FILLING: 0:56
HEATING: 25:11
STERILIZING: 30:00
VENTING: 2:53
DRYING: 26:23
TOTAL CYCLE: 01:25:23

06/07/2021 08:54 AM
ACCEPTED: X
REJECTED:
Operator: U01