
Midmark M11 - v1.0.0
Total Cycles: 592

Name: MIDMARK
Sterilizer ID: V2178198
01/11/2021 07:35 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 70.2 F 0.0
2:00 75.3 F 0.0
4:00 118.8 F 0.1
6:00 124.9 F 0.2
8:00 159.1 F 0.3
10:00 187.3 F 0.9
12:00 201.3 F 4.6
14:00 212.0 F 4.8
16:00 221.7 F 4.3
18:00 228.6 F 5.9
20:00 231.4 F 6.6
22:00 229.9 F 5.8
24:00 243.9 F 11.7

STERILIZING
mm:ss Degrees PSI
0:00 251.9 F 15.5
1:00 252.1 F 15.7
2:00 252.4 F 15.8
3:00 252.7 F 16.0
4:00 252.5 F 16.1
5:00 251.5 F 15.4
6:00 252.5 F 15.9
7:00 252.4 F 16.1
8:00 252.0 F 15.7
9:00 252.9 F 16.2
10:00 251.9 F 15.6
11:00 252.9 F 16.2
12:00 252.0 F 15.7
13:00 252.4 F 16.1
14:00 252.4 F 15.9
15:00 251.6 F 15.5
16:00 252.9 F 16.1
17:00 252.0 F 15.7
18:00 252.0 F 15.9
19:00 252.8 F 16.0
20:00 251.9 F 15.6
21:00 252.3 F 16.0
22:00 252.7 F 16.0
23:00 251.9 F 15.6
24:00 252.1 F 15.9
25:00 252.9 F 16.1
26:00 252.2 F 15.7
27:00 251.6 F 15.6
28:00 253.0 F 16.2
29:00 252.5 F 15.9
30:00 251.8 F 15.5

Min 251.4 F 15.4
Max 253.0 F 16.2

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 174.5 F 0.0
5:00 175.5 F 0.0
10:00 177.3 F 0.0
15:00 175.2 F 0.0
20:00 175.1 F 0.0
25:00 176.0 F 0.0
Dry Cycle Aborted

FILLING: 0:48
HEATING: 25:39
STERILIZING: 30:00
VENTING: 2:47
DRYING: 26:11
TOTAL CYCLE: 01:25:24

01/11/2021 09:00 AM
ACCEPTED: X
REJECTED:
Operator: U01