
Midmark M11 - v1.0.0
Total Cycles: 866

Name: RITTER 1
Sterilizer ID: V2178198
06/07/2021 11:25 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 191.8 F 0.0
2:00 186.6 F 4.1
4:00 201.7 F 3.4
6:00 217.9 F 5.7
8:00 221.5 F 4.0
10:00 226.6 F 5.0
12:00 237.3 F 9.4
14:00 250.4 F 15.7

STERILIZING
mm:ss Degrees PSI
0:00 251.5 F 16.2
1:00 252.6 F 16.6
2:00 251.6 F 16.3
3:00 252.7 F 16.7
4:00 251.8 F 16.2
5:00 252.3 F 16.8
6:00 252.6 F 16.7
7:00 252.2 F 16.4
8:00 251.6 F 16.2
9:00 252.5 F 16.8
10:00 252.7 F 16.8
11:00 252.5 F 16.6
12:00 252.2 F 16.4
13:00 251.8 F 16.2
14:00 251.8 F 16.4
15:00 252.5 F 16.8
16:00 252.8 F 16.8
17:00 252.7 F 16.7
18:00 252.5 F 16.6
19:00 252.3 F 16.5
20:00 252.1 F 16.4
21:00 251.9 F 16.2
22:00 251.7 F 16.2
23:00 252.2 F 16.7
24:00 252.6 F 16.9
25:00 252.8 F 16.8
26:00 252.7 F 16.8
27:00 252.6 F 16.7
28:00 252.4 F 16.5
29:00 252.3 F 16.4
30:00 252.0 F 16.3

Min 251.5 F 16.1
Max 253.0 F 16.9

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 199.0 F 0.0
5:00 225.4 F 0.0
10:00 212.8 F 0.0
15:00 195.6 F 0.0
Dry Cycle Aborted

FILLING: 1:41
HEATING: 14:46
STERILIZING: 30:00
VENTING: 2:56
DRYING: 15:41
TOTAL CYCLE: 01:05:04

06/07/2021 12:35 PM

ACCEPTED: X

REJECTED:

Operator: U01