
Midmark M11 - v1.0.0
Total Cycles: 590

Name: MIDMARK
Sterilizer ID: V2178198
12/30/2020 12:32 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 96.5 F 0.0
2:00 108.6 F 0.0
4:00 114.7 F 0.1
6:00 148.7 F 0.5
8:00 177.1 F 3.7
10:00 196.2 F 7.4
12:00 211.6 F 3.1
14:00 224.4 F 7.0
16:00 225.8 F 4.9
18:00 226.8 F 4.7
20:00 237.8 F 9.2
22:00 249.4 F 14.7

STERILIZING
mm:ss Degrees PSI
0:00 251.3 F 15.8
1:00 252.4 F 16.1
2:00 251.6 F 16.0
3:00 252.3 F 16.1
4:00 251.3 F 15.8
5:00 252.6 F 16.3
6:00 251.8 F 15.9
7:00 251.6 F 16.1
8:00 252.7 F 16.3
9:00 252.1 F 16.0
10:00 251.4 F 15.7
11:00 251.9 F 16.2
12:00 252.7 F 16.4
13:00 252.4 F 16.2
14:00 252.0 F 15.9
15:00 251.4 F 15.7
16:00 251.5 F 16.0
17:00 252.6 F 16.4
18:00 252.7 F 16.4
19:00 252.4 F 16.2
20:00 252.0 F 16.0
21:00 251.6 F 15.8
22:00 251.2 F 15.8
23:00 251.7 F 16.1
24:00 252.7 F 16.5
25:00 252.8 F 16.4
26:00 252.6 F 16.3
27:00 252.3 F 16.1
28:00 252.0 F 16.0
29:00 251.6 F 15.8
30:00 251.2 F 15.7

Min 251.2 F 15.7
Max 252.8 F 16.5

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 212.3 F 0.0
5:00 212.7 F 0.0
10:00 215.4 F 0.0
15:00 207.5 F 0.0
20:00 198.6 F 0.0
Dry Cycle Aborted

FILLING: 0:53

HEATING: 22:58
STERILIZING: 30:00
VENTING: 2:52
DRYING: 24:22
TOTAL CYCLE: 01:21:05

12/30/2020 01:53 PM
ACCEPTED: X
REJECTED:
Operator: U01