
Midmark M11 - v1.0.0
Total Cycles: 581

Name: MIDMARK
Sterilizer ID: V2178198
12/22/2020 09:58 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 105.6 F 0.0
2:00 107.5 F 0.0
4:00 130.6 F 0.1
6:00 161.1 F 1.5
8:00 188.0 F 5.1
10:00 208.1 F 2.4
12:00 218.4 F 5.2
14:00 222.0 F 3.9
16:00 231.4 F 6.8
18:00 231.9 F 6.9
20:00 244.4 F 12.4

STERILIZING
mm:ss Degrees PSI
0:00 251.5 F 16.0
1:00 251.8 F 16.0
2:00 252.2 F 16.2
3:00 252.6 F 16.5
4:00 252.0 F 16.4
5:00 251.9 F 16.0
6:00 252.5 F 16.4
7:00 251.8 F 16.2
8:00 252.2 F 16.2
9:00 252.5 F 16.6
10:00 252.1 F 16.2
11:00 252.5 F 16.6
12:00 252.2 F 16.2
13:00 252.1 F 16.4
14:00 252.4 F 16.3
15:00 251.6 F 16.0
16:00 252.6 F 16.5
17:00 251.9 F 16.1
18:00 252.6 F 16.6
19:00 252.3 F 16.3
20:00 251.6 F 16.0
21:00 252.8 F 16.6
22:00 252.1 F 16.2
23:00 251.8 F 16.2
24:00 252.7 F 16.5
25:00 252.0 F 16.1
26:00 252.0 F 16.4
27:00 252.6 F 16.5
28:00 251.9 F 16.1
29:00 252.2 F 16.5
30:00 252.6 F 16.4

Min 251.5 F 15.9
Max 252.9 F 16.6

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 178.3 F 0.0
5:00 202.6 F 0.0
10:00 182.8 F 0.0
15:00 180.9 F 0.0
20:00 201.4 F 0.0
25:00 219.4 F 0.0
30:00 226.7 F 0.0

Drying Complete

FILLING: 0:52
HEATING: 21:49
STERILIZING: 30:00
VENTING: 2:51
DRYING: 30:00
TOTAL CYCLE: 01:25:31
CYCLE COMPLETE

12/22/2020 11:23 AM
ACCEPTED: X
REJECTED:
Operator: U01