
Midmark M11 - v1.0.5
Total Cycles: 56

Name: RITTER2
Sterilizer ID: V2313228
05/31/2021 09:43 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 186.4 F 0.0
2:00 184.7 F 2.8
4:00 189.4 F 5.8
6:00 212.9 F 3.7
8:00 222.8 F 5.2
10:00 231.8 F 7.0
12:00 234.5 F 8.0
14:00 247.9 F 14.2

STERILIZING
mm:ss Degrees PSI
0:00 251.4 F 16.0
1:00 252.2 F 16.2
2:00 252.3 F 16.6
3:00 251.9 F 16.1
4:00 252.3 F 16.6
5:00 252.3 F 16.3
6:00 251.5 F 16.0
7:00 252.7 F 16.7
8:00 252.3 F 16.4
9:00 251.7 F 16.1
10:00 251.6 F 16.2
11:00 252.8 F 16.7
12:00 252.5 F 16.5
13:00 252.2 F 16.3
14:00 251.8 F 16.1
15:00 251.4 F 16.0
16:00 251.8 F 16.3
17:00 252.6 F 16.7
18:00 252.8 F 16.6
19:00 252.6 F 16.5
20:00 252.4 F 16.4
21:00 252.1 F 16.3
22:00 252.0 F 16.2
23:00 251.8 F 16.1
24:00 251.5 F 16.0
25:00 251.4 F 16.0
26:00 251.8 F 16.4
27:00 252.4 F 16.7
28:00 252.8 F 16.7
29:00 252.8 F 16.6
30:00 252.7 F 16.6

Min 251.4 F 15.9
Max 252.9 F 16.7

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 185.7 F 0.0
5:00 202.5 F 0.0
10:00 201.6 F 0.0
15:00 199.8 F 0.0
20:00 201.3 F 0.0
25:00 204.7 F 0.0
Dry Cycle Aborted

FILLING: 1:57
HEATING: 15:09
STERILIZING: 30:00
VENTING: 2:51

DRYING: 27:28
TOTAL CYCLE: 01:17:25

05/31/2021 11:00 AM
ACCEPTED: X
REJECTED:
Operator: U01