

-----  
Midmark M11 - v1.0.5  
Total Cycles: 210

Name: RITTER2  
Sterilizer ID: V2313228  
07/14/2021 02:45 PM

BEGIN PACKS CYCLE  
Temp: 250 Degrees F  
Time: 30 Minutes  
Dry: 30 Minutes

FILLING CHAMBER

HEATING  
mm:ss Degrees PSI  
0:00 154.8 F 0.0  
2:00 163.4 F 2.5  
4:00 191.1 F 6.4  
6:00 209.5 F 2.3  
8:00 216.0 F 6.8  
10:00 221.8 F 4.8  
12:00 225.8 F 4.9  
14:00 236.7 F 9.5  
16:00 248.6 F 15.1

STERILIZING  
mm:ss Degrees PSI  
0:00 251.4 F 16.7  
1:00 252.5 F 17.1  
2:00 251.3 F 16.7  
3:00 252.6 F 17.2  
4:00 251.7 F 16.8  
5:00 251.4 F 17.0  
6:00 252.6 F 17.2  
7:00 252.1 F 17.0  
8:00 251.6 F 16.8  
9:00 251.2 F 16.7  
10:00 251.5 F 17.1  
11:00 252.0 F 17.3  
12:00 252.7 F 17.3  
13:00 252.4 F 17.2  
14:00 252.5 F 17.2  
15:00 252.3 F 17.1  
16:00 252.2 F 17.1  
17:00 252.1 F 17.0  
18:00 252.1 F 17.0  
19:00 252.0 F 17.0  
20:00 252.0 F 17.0  
21:00 251.9 F 17.0  
22:00 252.0 F 17.0  
23:00 251.9 F 17.0  
24:00 252.0 F 17.0  
25:00 252.0 F 17.0  
26:00 252.0 F 17.0  
27:00 252.0 F 17.0  
28:00 252.0 F 17.0  
29:00 252.0 F 17.0  
30:00 252.0 F 17.0

Min 251.0 F 16.6  
Max 252.9 F 17.4

VENTING CHAMBER

DRYING  
Drying Start  
mm:ss Degrees PSI  
0:00 201.4 F 0.0  
Dry Cycle Aborted

FILLING: 1:40  
HEATING: 17:08  
STERILIZING: 30:00  
VENTING: 3:05  
DRYING: 3:57  
TOTAL CYCLE: 00:55:49

07/14/2021 03:40 PM

ACCEPTED: X  
REJECTED:  
Operator: U01