
Midmark M11 - v1.0.0
Total Cycles: 576

Name: MIDMARK
Sterilizer ID: V2178198
12/21/2020 09:25 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 156.3 F 0.0
2:00 156.3 F 1.3
4:00 163.9 F 4.3
6:00 191.2 F 1.6
8:00 214.5 F 4.2
10:00 222.5 F 4.4
12:00 231.2 F 6.9
14:00 232.4 F 7.1
16:00 245.7 F 13.0

STERILIZING
mm:ss Degrees PSI
0:00 251.6 F 16.0
1:00 252.4 F 16.4
2:00 251.7 F 16.1
3:00 252.4 F 16.4
4:00 251.6 F 15.9
5:00 252.7 F 16.6
6:00 252.2 F 16.3
7:00 251.7 F 15.9
8:00 252.2 F 16.5
9:00 252.6 F 16.6
10:00 252.2 F 16.3
11:00 251.8 F 16.0
12:00 251.6 F 16.0
13:00 252.5 F 16.6
14:00 252.7 F 16.6
15:00 252.4 F 16.5
16:00 252.0 F 16.3
17:00 251.8 F 16.1
18:00 251.5 F 15.9
19:00 251.8 F 16.2
20:00 252.7 F 16.6
21:00 252.7 F 16.6
22:00 252.5 F 16.5
23:00 252.3 F 16.4
24:00 252.0 F 16.2
25:00 251.8 F 16.1
26:00 251.6 F 15.9
27:00 251.5 F 16.0
28:00 252.0 F 16.4
29:00 252.7 F 16.7
30:00 252.8 F 16.7

Min 251.4 F 15.9
Max 252.9 F 16.7

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 184.1 F 0.0
5:00 184.1 F 0.0
10:00 179.5 F 0.0
15:00 179.5 F 0.0
20:00 177.7 F 0.0
25:00 176.7 F 0.0
Dry Cycle Aborted

FILLING: 0:59
HEATING: 17:35
STERILIZING: 30:00

VENTING: 2:53
DRYING: 26:39
TOTAL CYCLE: 01:18:05

12/21/2020 10:43 AM

ACCEPTED:

REJECTED:

Reason: -----

Operator: -----