
Midmark M11 - v1.0.0
Total Cycles: 588

Name: MIDMARK
Sterilizer ID: V2178198
12/30/2020 07:39 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 73.7 F 0.0
2:00 81.0 F 0.0
4:00 91.2 F 0.1
6:00 129.6 F 0.1
8:00 162.4 F 0.3
10:00 187.6 F 1.3
12:00 211.1 F 4.7
14:00 216.2 F 5.4
16:00 222.3 F 4.4
18:00 227.4 F 5.5
20:00 231.7 F 6.7
22:00 229.9 F 5.8
24:00 243.6 F 11.6

STERILIZING
mm:ss Degrees PSI
0:00 251.3 F 15.6
1:00 252.4 F 15.9
2:00 252.1 F 16.0
3:00 251.9 F 15.7
4:00 252.5 F 16.2
5:00 252.0 F 15.7
6:00 251.9 F 15.9
7:00 252.6 F 16.0
8:00 251.6 F 15.5
9:00 252.3 F 16.2
10:00 252.5 F 16.0
11:00 251.8 F 15.6
12:00 251.9 F 15.9
13:00 252.7 F 16.1
14:00 252.2 F 15.8
15:00 251.5 F 15.5
16:00 252.2 F 16.1
17:00 252.7 F 16.2
18:00 252.3 F 15.9
19:00 251.7 F 15.6
20:00 251.6 F 15.8
21:00 252.7 F 16.2
22:00 252.6 F 16.1
23:00 252.1 F 15.8
24:00 251.6 F 15.5
25:00 251.6 F 15.8
26:00 252.7 F 16.3
27:00 252.6 F 16.1
28:00 252.2 F 15.8
29:00 251.7 F 15.6
30:00 251.6 F 15.8

Min 251.3 F 15.4
Max 252.8 F 16.3

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 178.1 F 0.0
5:00 220.6 F 0.0
10:00 216.7 F 0.0
15:00 223.0 F 0.0
20:00 187.8 F 0.0
25:00 184.2 F 0.0
Dry Cycle Aborted

FILLING: 0:59
HEATING: 25:55
STERILIZING: 30:00
VENTING: 2:49
DRYING: 26:20
TOTAL CYCLE: 01:26:03

12/30/2020 09:20 AM
ACCEPTED: X
REJECTED:
Operator: U01