
Midmark M11 - v1.0.5
Total Cycles: 4

Name: RITTER2
Sterilizer ID: V2313228
04/21/2021 12:10 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 148.4 F 0.0
2:00 148.9 F 2.9
4:00 175.2 F 7.5
6:00 203.9 F 3.1
8:00 217.8 F 5.6
10:00 229.5 F 7.0
12:00 232.0 F 7.2
14:00 245.9 F 13.4

STERILIZING
mm:ss Degrees PSI
0:00 252.0 F 16.3
1:00 252.3 F 16.5
2:00 252.8 F 16.7
3:00 252.6 F 17.0
4:00 251.6 F 16.3
5:00 252.5 F 16.6
6:00 252.2 F 16.8
7:00 252.2 F 16.5
8:00 252.6 F 17.0
9:00 252.2 F 16.5
10:00 252.5 F 17.0
11:00 252.5 F 16.6
12:00 251.8 F 16.6
13:00 252.7 F 16.8
14:00 251.8 F 16.3
15:00 252.9 F 17.0
16:00 252.4 F 16.6
17:00 251.7 F 16.2
18:00 252.9 F 17.0
19:00 252.4 F 16.6
20:00 251.6 F 16.2
21:00 252.9 F 17.0
22:00 252.4 F 16.6
23:00 251.7 F 16.2
24:00 252.8 F 17.0
25:00 252.4 F 16.6
26:00 251.7 F 16.2
27:00 252.8 F 17.0
28:00 252.4 F 16.6
29:00 251.7 F 16.3
30:00 252.7 F 17.0

Min 251.6 F 16.2
Max 253.1 F 17.0

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 171.7 F 0.0
5:00 167.7 F 0.0
10:00 161.4 F 0.0
15:00 177.5 F 0.0
20:00 161.5 F 0.0
25:00 160.0 F 0.0
Dry Cycle Aborted

FILLING: 1:22
HEATING: 15:20
STERILIZING: 30:00
VENTING: 2:45

DRYING: 25:42
TOTAL CYCLE: 01:15:08

04/21/2021 01:25 PM
ACCEPTED: X
REJECTED:
Operator: U01