
Midmark M11 - v1.0.0
Total Cycles: 2079

Name: RITTER 1
Sterilizer ID: V2178198
07/04/2022 07:30 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 72.3 F 0.0
2:00 76.4 F 0.0
4:00 90.9 F 0.0
6:00 123.6 F 0.8
8:00 145.2 F 1.7
10:00 162.9 F 3.4
12:00 177.9 F 5.4
14:00 200.5 F 3.0
16:00 216.7 F 4.1
18:00 225.3 F 6.7
20:00 226.3 F 5.4
22:00 224.7 F 4.1
24:00 237.4 F 9.1
26:00 248.7 F 14.5

STERILIZING
mm:ss Degrees PSI
0:00 251.8 F 16.0
1:00 252.5 F 16.5
2:00 251.2 F 15.9
3:00 252.4 F 16.4
4:00 251.2 F 16.0
5:00 252.4 F 16.5
6:00 251.4 F 15.9
7:00 252.3 F 16.6
8:00 251.9 F 16.2
9:00 251.3 F 16.1
10:00 252.4 F 16.5
11:00 251.7 F 16.1
12:00 251.3 F 16.1
13:00 252.4 F 16.5
14:00 251.8 F 16.2
15:00 251.2 F 16.0
16:00 252.5 F 16.6
17:00 252.1 F 16.3
18:00 251.4 F 15.9
19:00 252.0 F 16.5
20:00 252.3 F 16.4
21:00 251.7 F 16.1
22:00 251.3 F 16.1
23:00 252.6 F 16.6
24:00 252.1 F 16.3
25:00 251.4 F 16.0
26:00 251.5 F 16.3
27:00 252.5 F 16.5
28:00 252.1 F 16.3
29:00 251.5 F 16.0
30:00 251.7 F 16.4

Min 251.1 F 15.9
Max 252.7 F 16.6

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 185.7 F 0.0
5:00 191.1 F 0.0
10:00 185.5 F 0.0
15:00 195.2 F 0.0
20:00 207.7 F 0.0
25:00 210.4 F 0.0

30:00 213.6 F 0.0
Drying Complete

FILLING: 0:45
HEATING: 26:60
STERILIZING: 30:00
VENTING: 2:52
DRYING: 30:00
TOTAL CYCLE: 01:30:37
CYCLE COMPLETE

07/04/2022 09:00 AM

ACCEPTED:

REJECTED:

Reason: _____
Operator: _____