
Midmark M11 - v1.0.0
Total Cycles: 595

Name: MIDMARK
Sterilizer ID: V2178198
01/11/2021 02:57 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 95.8 F 0.0
2:00 96.2 F 0.0
4:00 116.4 F 0.4
6:00 149.5 F 2.6
8:00 189.9 F 6.2
10:00 203.5 F 2.1
12:00 217.9 F 6.2
14:00 220.8 F 3.6
16:00 229.0 F 5.8
18:00 234.0 F 7.7
20:00 246.5 F 13.3

STERILIZING
mm:ss Degrees PSI
0:00 251.3 F 15.9
1:00 252.4 F 16.2
2:00 252.0 F 16.3
3:00 252.0 F 16.0
4:00 252.4 F 16.5
5:00 252.2 F 16.1
6:00 251.6 F 16.1
7:00 252.6 F 16.4
8:00 251.6 F 15.9
9:00 252.1 F 16.4
10:00 252.6 F 16.3
11:00 251.9 F 16.0
12:00 251.5 F 16.1
13:00 252.7 F 16.5
14:00 252.4 F 16.2
15:00 251.7 F 15.9
16:00 251.6 F 16.1
17:00 252.8 F 16.5
18:00 252.5 F 16.3
19:00 252.0 F 16.0
20:00 251.4 F 15.8
21:00 252.1 F 16.4
22:00 252.8 F 16.5
23:00 252.4 F 16.3
24:00 252.0 F 16.0
25:00 251.3 F 15.9
26:00 252.1 F 16.5
27:00 252.8 F 16.5
28:00 252.5 F 16.3
29:00 252.1 F 16.1
30:00 251.6 F 15.9

Min 251.3 F 15.8
Max 252.9 F 16.6

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 187.8 F 0.0
5:00 202.1 F 0.0
10:00 186.4 F 0.0
Dry Cycle Aborted

FILLING: 0:48
HEATING: 21:29
STERILIZING: 30:00
VENTING: 2:49

DRYING: 12:01
TOTAL CYCLE: 01:07:08

01/11/2021 04:04 PM
ACCEPTED: X
REJECTED:
Operator: U01