
Midmark M11 - v1.0.0
Total Cycles: 1033

Name: RITTER 1
Sterilizer ID: V2178198
07/20/2021 12:06 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 99.6 F 0.0
2:00 98.8 F 0.0
4:00 103.8 F 0.0
6:00 131.1 F 0.1
8:00 163.5 F 0.9
10:00 184.2 F 3.4
12:00 202.5 F 6.1
14:00 213.3 F 2.3
16:00 227.7 F 6.1
18:00 223.1 F 4.0
20:00 227.7 F 5.6
22:00 234.9 F 7.9
24:00 245.9 F 12.9
26:00 252.2 F 16.3

STERILIZING
mm:ss Degrees PSI
0:00 251.8 F 16.2
1:00 251.8 F 16.4
2:00 251.9 F 16.5
3:00 251.9 F 16.5
4:00 251.8 F 16.5
5:00 251.5 F 16.3
6:00 252.0 F 16.1
7:00 252.6 F 16.3
8:00 252.0 F 16.5
9:00 251.6 F 16.4
10:00 251.8 F 16.2
11:00 252.1 F 16.4
12:00 251.7 F 16.5
13:00 251.8 F 16.2
14:00 251.9 F 16.4
15:00 251.4 F 16.3
16:00 252.0 F 16.3
17:00 251.8 F 16.5
18:00 251.3 F 16.1
19:00 251.6 F 16.4
20:00 251.5 F 16.3
21:00 251.7 F 16.3
22:00 251.7 F 16.5
23:00 251.7 F 16.2
24:00 251.8 F 16.5
25:00 251.3 F 16.2
26:00 251.9 F 16.4
27:00 251.3 F 16.4
28:00 251.3 F 16.3
29:00 251.8 F 16.5
30:00 251.6 F 16.2

Min 251.1 F 16.1
Max 252.7 F 16.5

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 195.5 F 0.0
5:00 182.1 F 0.0
10:00 197.3 F 0.0
15:00 199.8 F 0.0
20:00 207.0 F 0.0
25:00 179.7 F 0.0

30:00 173.6 F 0.0
Drying Complete

FILLING: 0:54
HEATING: 26:23
STERILIZING: 30:00
VENTING: 3:07
DRYING: 30:00
TOTAL CYCLE: 01:30:24
CYCLE COMPLETE

07/20/2021 01:43 PM
ACCEPTED: X
REJECTED:
Operator: U01