
Midmark M11 - v1.0.0
Total Cycles: 655

Name: MIDMARK
Sterilizer ID: V2178198
02/09/2021 07:42 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

mm:ss	Degrees	PSI
0:00	71.1 F	0.0
2:00	72.1 F	0.0
4:00	73.9 F	0.1
6:00	102.0 F	0.2
8:00	141.8 F	0.2
10:00	183.2 F	1.5
12:00	217.1 F	4.6
14:00	231.2 F	6.8
16:00	225.0 F	4.2
18:00	226.9 F	5.8
20:00	230.4 F	6.3
22:00	227.9 F	5.1
24:00	241.7 F	10.8
26:00	252.6 F	16.0

STERILIZING

mm:ss	Degrees	PSI
0:00	251.5 F	15.6
1:00	252.8 F	16.0
2:00	251.8 F	15.8
3:00	252.3 F	16.0
4:00	251.6 F	15.8
5:00	252.9 F	16.1
6:00	251.7 F	15.6
7:00	252.6 F	16.3
8:00	252.7 F	16.0
9:00	251.6 F	15.6
10:00	252.2 F	16.2
11:00	252.7 F	16.2
12:00	251.9 F	15.8
13:00	251.5 F	15.6
14:00	252.7 F	16.3
15:00	252.7 F	16.2
16:00	252.0 F	15.9
17:00	251.5 F	15.6
18:00	251.8 F	15.9
19:00	252.9 F	16.3
20:00	252.6 F	16.2
21:00	252.3 F	15.9
22:00	251.5 F	15.6
23:00	251.5 F	15.7
24:00	252.7 F	16.3
25:00	252.8 F	16.2
26:00	252.3 F	16.0
27:00	252.2 F	15.8
28:00	251.4 F	15.6
29:00	251.7 F	15.9
30:00	252.8 F	16.3

Min 251.3 F 15.5
Max 253.1 F 16.3

VENTING CHAMBER

DRYING
Drying Start

mm:ss	Degrees	PSI
0:00	192.4 F	0.0
5:00	214.7 F	0.0
10:00	192.7 F	0.0
15:00	197.2 F	0.0
20:00	177.8 F	0.0

Dry Cycle Aborted

FILLING: 0:49
HEATING: 26:17
STERILIZING: 30:00
VENTING: 2:48
DRYING: 20:07
TOTAL CYCLE: 01:20:01

02/09/2021 09:02 AM
ACCEPTED: X
REJECTED:
Operator: U01