
Midmark M11 - v1.0.0
Total Cycles: 591

Name: MIDMARK
Sterilizer ID: V2178198
12/30/2020 02:39 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 125.3 F 0.0
2:00 124.2 F 0.6
4:00 141.6 F 2.1
6:00 180.0 F 5.4
8:00 203.4 F 1.7
10:00 214.3 F 6.0
12:00 219.9 F 3.5
14:00 228.4 F 5.7
16:00 233.5 F 7.6
18:00 245.9 F 13.2

STERILIZING
mm:ss Degrees PSI
0:00 251.3 F 16.0
1:00 252.0 F 16.1
2:00 252.4 F 16.4
3:00 251.7 F 16.4
4:00 251.8 F 16.1
5:00 252.5 F 16.6
6:00 251.6 F 16.0
7:00 252.6 F 16.6
8:00 251.8 F 16.1
9:00 252.0 F 16.5
10:00 252.4 F 16.3
11:00 251.2 F 16.0
12:00 252.6 F 16.6
13:00 252.1 F 16.3
14:00 251.2 F 16.0
15:00 252.6 F 16.6
16:00 252.3 F 16.3
17:00 251.3 F 16.0
18:00 252.4 F 16.6
19:00 252.6 F 16.5
20:00 251.8 F 16.2
21:00 251.3 F 16.1
22:00 252.6 F 16.6
23:00 252.4 F 16.4
24:00 251.6 F 16.1
25:00 251.4 F 16.3
26:00 252.7 F 16.6
27:00 252.5 F 16.4
28:00 251.7 F 16.1
29:00 251.5 F 16.2
30:00 252.7 F 16.6

Min 251.2 F 15.9
Max 252.8 F 16.6

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 204.9 F 0.0
5:00 215.6 F 0.0
10:00 187.7 F 0.0
15:00 182.2 F 0.0
20:00 181.3 F 0.0
25:00 197.4 F 0.0
30:00 226.6 F 0.0
Drying Complete

FILLING: 1:02

HEATING: 19:31
STERILIZING: 30:00
VENTING: 2:53
DRYING: 30:00
TOTAL CYCLE: 01:23:26
CYCLE COMPLETE

12/30/2020 04:02 PM
ACCEPTED: X
REJECTED:
Operator: U01