
Midmark M11 - v1.0.5
Total Cycles: 14

Name: RITTER2
Sterilizer ID: V2313228
04/26/2021 06:19 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 71.7 F 0.0
2:00 72.2 F 0.2
4:00 77.1 F 1.1
6:00 100.9 F 3.6
8:00 131.3 F 6.8
10:00 174.0 F 0.8
12:00 193.6 F 4.1
14:00 206.1 F 5.7
16:00 219.9 F 4.4
18:00 226.5 F 5.7
20:00 232.5 F 7.6
22:00 243.9 F 12.7

STERILIZING
mm:ss Degrees PSI
0:00 252.1 F 16.5
1:00 252.4 F 16.8
2:00 252.8 F 17.0
3:00 252.1 F 16.9
4:00 251.8 F 16.4
5:00 252.4 F 16.9
6:00 252.2 F 17.0
7:00 252.1 F 16.6
8:00 252.5 F 17.1
9:00 252.0 F 16.6
10:00 252.6 F 17.1
11:00 252.1 F 16.6
12:00 252.4 F 17.0
13:00 252.4 F 16.8
14:00 251.8 F 16.5
15:00 252.7 F 17.1
16:00 252.1 F 16.7
17:00 251.8 F 16.7
18:00 252.7 F 17.1
19:00 252.3 F 16.7
20:00 251.8 F 16.5
21:00 252.7 F 17.2
22:00 252.5 F 16.9
23:00 252.0 F 16.6
24:00 252.2 F 16.9
25:00 252.7 F 17.1
26:00 252.3 F 16.8
27:00 252.0 F 16.5
28:00 252.4 F 17.0
29:00 252.7 F 17.1
30:00 252.3 F 16.8

Min 251.6 F 16.4
Max 252.9 F 17.2

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 173.6 F 0.0
5:00 195.6 F 0.0
10:00 202.1 F 0.0
15:00 180.2 F 0.0
20:00 175.3 F 0.0
25:00 161.1 F 0.0
Dry Cycle Aborted

FILLING: 0:57
HEATING: 23:48
STERILIZING: 30:00
VENTING: 2:48
DRYING: 26:37
TOTAL CYCLE: 01:24:10

04/26/2021 07:43 AM
ACCEPTED: X
REJECTED:
Operator: U01