
Midmark M11 - v1.0.0
Total Cycles: 862

Name: RITTER 1
Sterilizer ID: V2178198
06/07/2021 06:37 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 70.1 F 0.0
2:00 76.0 F 0.0
4:00 106.1 F 0.1
6:00 127.1 F 0.1
8:00 156.2 F 0.2
10:00 177.3 F 0.2
12:00 198.9 F 0.7
14:00 213.5 F 3.9
16:00 222.0 F 6.4
18:00 218.3 F 3.0
20:00 230.6 F 7.1
22:00 228.7 F 5.6
24:00 226.2 F 4.4
26:00 240.2 F 10.2
28:00 251.7 F 15.8

STERILIZING
mm:ss Degrees PSI
0:00 251.5 F 15.7
1:00 252.6 F 16.1
2:00 251.5 F 15.7
3:00 252.7 F 16.2
4:00 251.8 F 15.7
5:00 252.2 F 16.2
6:00 252.7 F 16.2
7:00 251.9 F 15.8
8:00 251.6 F 15.8
9:00 252.8 F 16.4
10:00 252.7 F 16.2
11:00 252.2 F 15.9
12:00 251.7 F 15.7
13:00 251.6 F 15.8
14:00 252.6 F 16.3
15:00 252.9 F 16.3
16:00 252.6 F 16.2
17:00 252.3 F 16.0
18:00 251.9 F 15.8
19:00 251.5 F 15.6
20:00 251.8 F 15.9
21:00 252.7 F 16.4
22:00 252.9 F 16.4
23:00 252.7 F 16.2
24:00 252.5 F 16.1
25:00 252.1 F 15.9
26:00 251.8 F 15.7
27:00 251.5 F 15.7
28:00 251.9 F 16.0
29:00 252.6 F 16.4
30:00 252.9 F 16.4

Min 251.5 F 15.6
Max 252.9 F 16.4

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 179.5 F 0.0
5:00 205.2 F 0.0
10:00 174.6 F 0.0
15:00 159.3 F 0.0
20:00 178.0 F 0.0

25:00 181.4 F 0.0
Dry Cycle Aborted

FILLING: 0:51
HEATING: 28:36
STERILIZING: 30:00
VENTING: 2:59
DRYING: 26:43
TOTAL CYCLE: 01:29:09

06/07/2021 08:06 AM
ACCEPTED: X
REJECTED:
Operator: U01