

-----  
Midmark M11 - v1.0.5  
Total Cycles: 1774

Name: RITTER2  
Sterilizer ID: V2313228  
11/28/2022 09:02 AM

BEGIN PACKS CYCLE  
Temp: 250 Degrees F  
Time: 30 Minutes  
Dry: 30 Minutes

FILLING CHAMBER

HEATING  
mm:ss Degrees PSI  
0:00 187.2 F 0.0  
2:00 182.8 F 0.8  
4:00 186.9 F 4.9  
6:00 212.0 F 3.0  
8:00 224.2 F 5.6  
10:00 229.6 F 6.8  
12:00 232.9 F 7.9  
14:00 247.2 F 14.3

STERILIZING  
mm:ss Degrees PSI  
0:00 251.2 F 16.6  
1:00 252.6 F 17.1  
2:00 251.3 F 16.6  
3:00 252.9 F 17.2  
4:00 252.1 F 16.8  
5:00 251.3 F 16.7  
6:00 252.5 F 17.3  
7:00 252.6 F 17.1  
8:00 252.2 F 16.8  
9:00 251.5 F 16.6  
10:00 251.2 F 16.7  
11:00 252.1 F 17.2  
12:00 252.5 F 17.3  
13:00 252.9 F 17.2  
14:00 252.5 F 17.1  
15:00 252.5 F 17.0  
16:00 252.3 F 17.0  
17:00 252.2 F 16.9  
18:00 251.9 F 16.8  
19:00 251.7 F 16.7  
20:00 251.5 F 16.6  
21:00 251.3 F 16.6  
22:00 251.3 F 16.6  
23:00 251.2 F 16.6  
24:00 251.3 F 16.8  
25:00 251.5 F 16.9  
26:00 251.8 F 17.1  
27:00 251.9 F 17.2  
28:00 252.3 F 17.3  
29:00 252.5 F 17.3  
30:00 252.9 F 17.3

Min 251.2 F 16.5  
Max 253.1 F 17.4

VENTING CHAMBER

DRYING  
Drying Start  
mm:ss Degrees PSI  
0:00 197.6 F 0.0  
5:00 206.1 F 0.0  
10:00 211.5 F 0.0  
15:00 214.5 F 0.0  
20:00 218.7 F 0.0  
25:00 220.0 F 0.0  
30:00 220.3 F 0.0  
Drying Complete

FILLING: 1:40  
HEATING: 15:19  
STERILIZING: 30:00

VENTING: 2:47  
DRYING: 30:00  
TOTAL CYCLE: 01:19:47  
CYCLE COMPLETE

11/28/2022 10:21 AM  
ACCEPTED:  
REJECTED:  
Reason: \_\_\_\_\_  
Operator: \_\_\_\_\_