

-----  
Midmark M11 - v1.0.0  
Total Cycles: 575

Name: MIDMARK  
Sterilizer ID: V2178198  
12/21/2020 07:38 AM

BEGIN PACKS CYCLE  
Temp: 250 Degrees F  
Time: 30 Minutes  
Dry: 30 Minutes

FILLING CHAMBER

HEATING  
mm:ss Degrees PSI  
0:00 71.2 F 0.0  
2:00 75.1 F 0.0  
4:00 101.8 F 0.0  
6:00 143.9 F 0.2  
8:00 149.6 F 0.3  
10:00 177.6 F 1.8  
12:00 198.4 F 4.8  
14:00 208.8 F 3.8  
16:00 217.3 F 4.3  
18:00 226.0 F 6.1  
20:00 227.5 F 5.6  
22:00 225.4 F 4.1  
24:00 238.0 F 9.4  
26:00 249.4 F 14.8

STERILIZING  
mm:ss Degrees PSI  
0:00 251.4 F 15.9  
1:00 252.2 F 16.1  
2:00 252.2 F 16.4  
3:00 251.8 F 15.9  
4:00 252.6 F 16.5  
5:00 251.8 F 15.9  
6:00 252.5 F 16.5  
7:00 252.2 F 16.1  
8:00 251.5 F 16.0  
9:00 252.7 F 16.5  
10:00 252.0 F 16.0  
11:00 251.5 F 16.0  
12:00 252.7 F 16.5  
13:00 252.2 F 16.2  
14:00 251.5 F 15.8  
15:00 252.5 F 16.5  
16:00 252.7 F 16.4  
17:00 252.1 F 16.1  
18:00 251.5 F 15.8  
19:00 252.2 F 16.4  
20:00 252.7 F 16.5  
21:00 252.3 F 16.2  
22:00 251.8 F 15.9  
23:00 251.6 F 16.1  
24:00 252.8 F 16.5  
25:00 252.6 F 16.4  
26:00 252.2 F 16.2  
27:00 251.7 F 15.9  
28:00 251.6 F 16.1  
29:00 252.8 F 16.6  
30:00 252.7 F 16.4

Min 251.4 F 15.8  
Max 252.9 F 16.6

VENTING CHAMBER

DRYING  
Drying Start  
mm:ss Degrees PSI  
0:00 207.9 F 0.0  
5:00 193.4 F 0.0  
10:00 188.1 F 0.0  
15:00 186.1 F 0.0  
20:00 204.2 F 0.0  
25:00 208.7 F 0.0

Dry Cycle Aborted

FILLING: 0:49  
HEATING: 26:59  
STERILIZING: 30:00  
VENTING: 2:52  
DRYING: 29:08  
TOTAL CYCLE: 01:29:48

12/21/2020 09:07 AM

ACCEPTED:

REJECTED:

Reason: \_\_\_\_\_  
Operator: \_\_\_\_\_