
Midmark M11 - v1.0.5
Total Cycles: 58

Name: RITTER2
Sterilizer ID: V2313228
05/31/2021 12:28 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 130.9 F 0.0
2:00 134.0 F 1.6
4:00 147.1 F 5.0
6:00 187.6 F 1.1
8:00 212.5 F 5.0
10:00 217.3 F 3.8
12:00 228.7 F 6.6
14:00 232.2 F 7.3
16:00 245.1 F 13.1

STERILIZING
mm:ss Degrees PSI
0:00 251.7 F 16.3
1:00 252.6 F 16.8
2:00 251.0 F 16.3
3:00 253.0 F 17.1
4:00 253.1 F 17.1
5:00 252.5 F 17.2
6:00 252.0 F 16.8
7:00 252.6 F 16.9
8:00 252.6 F 17.2
9:00 252.3 F 16.7
10:00 253.1 F 17.1
11:00 251.8 F 16.8
12:00 252.6 F 17.1
13:00 251.8 F 16.7
14:00 253.1 F 17.1
15:00 252.3 F 16.8
16:00 253.1 F 17.2
17:00 252.5 F 16.9
18:00 252.1 F 17.0
19:00 253.3 F 17.1
20:00 252.4 F 16.8
21:00 252.2 F 17.2
22:00 252.6 F 17.1
23:00 252.2 F 16.7
24:00 252.5 F 17.3
25:00 252.8 F 17.0
26:00 252.2 F 16.7
27:00 252.6 F 17.2
28:00 253.0 F 17.0
29:00 252.2 F 16.7
30:00 252.5 F 17.3

Min 250.9 F 16.2
Max 253.3 F 17.3

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 195.5 F 0.0
5:00 173.0 F 0.0
10:00 162.6 F 0.0
15:00 157.2 F 0.0
20:00 167.4 F 0.0
25:00 161.5 F 0.0
Dry Cycle Aborted

FILLING: 1:16
HEATING: 17:32
STERILIZING: 30:00

VENTING: 2:50
DRYING: 26:10
TOTAL CYCLE: 01:17:48

05/31/2021 01:45 PM
ACCEPTED: X
REJECTED:
Operator: U01