
Midmark M11 - v1.0.5
Total Cycles: 73

Name: RITTER2
Sterilizer ID: V2313228
06/07/2021 01:53 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 192.3 F 0.0
2:00 192.3 F 4.7
4:00 212.5 F 3.2
6:00 224.9 F 5.8
8:00 229.8 F 6.5
10:00 236.6 F 9.2
12:00 250.9 F 16.0

STERILIZING
mm:ss Degrees PSI
0:00 251.8 F 16.4
1:00 251.9 F 16.3
2:00 252.2 F 16.5
3:00 252.5 F 16.8
4:00 252.3 F 16.9
5:00 252.1 F 16.5
6:00 252.7 F 17.0
7:00 252.0 F 16.5
8:00 252.8 F 17.0
9:00 252.1 F 16.6
10:00 252.2 F 16.9
11:00 252.6 F 16.8
12:00 252.0 F 16.5
13:00 252.2 F 17.0
14:00 252.6 F 16.9
15:00 252.2 F 16.6
16:00 251.7 F 16.5
17:00 252.7 F 17.1
18:00 252.5 F 16.9
19:00 252.1 F 16.6
20:00 251.7 F 16.4
21:00 252.7 F 17.1
22:00 252.6 F 16.9
23:00 252.3 F 16.7
24:00 251.9 F 16.5
25:00 251.8 F 16.5
26:00 252.6 F 17.1
27:00 252.7 F 17.0
28:00 252.4 F 16.8
29:00 252.2 F 16.6
30:00 251.8 F 16.4

Min 251.6 F 16.3
Max 252.8 F 17.1

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 194.6 F 0.0
5:00 217.2 F 0.0
10:00 220.0 F 0.0
15:00 216.0 F 0.0
Dry Cycle Aborted

FILLING: 2:53
HEATING: 12:35
STERILIZING: 30:00
VENTING: 2:50
DRYING: 18:60
TOTAL CYCLE: 01:07:18

06/07/2021 03:00 PM
ACCEPTED: X
REJECTED:
Operator: U01