
Midmark M11 - v1.0.0
Total Cycles: 690

Name: MIDMARK
Sterilizer ID: V2178198
02/23/2021 09:21 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 111.5 F 0.0
2:00 113.3 F 0.5
4:00 123.7 F 2.1
6:00 158.9 F 5.1
8:00 192.6 F 1.0
10:00 209.6 F 3.8
12:00 221.9 F 6.0
14:00 227.0 F 5.8
16:00 224.8 F 4.5
18:00 239.0 F 10.1
20:00 250.6 F 15.7

STERILIZING
mm:ss Degrees PSI
0:00 251.3 F 16.1
1:00 252.3 F 16.5
2:00 251.6 F 16.4
3:00 252.3 F 16.4
4:00 251.4 F 16.3
5:00 252.5 F 16.6
6:00 251.7 F 16.2
7:00 251.9 F 16.6
8:00 252.4 F 16.6
9:00 251.8 F 16.2
10:00 251.3 F 16.3
11:00 252.6 F 16.8
12:00 252.4 F 16.5
13:00 251.8 F 16.2
14:00 251.2 F 16.2
15:00 252.2 F 16.8
16:00 252.5 F 16.7
17:00 252.2 F 16.5
18:00 251.9 F 16.3
19:00 251.3 F 16.1
20:00 251.5 F 16.4
21:00 252.5 F 16.8
22:00 252.5 F 16.7
23:00 252.3 F 16.6
24:00 251.9 F 16.4
25:00 251.7 F 16.2
26:00 251.3 F 16.1
27:00 251.7 F 16.5
28:00 252.6 F 16.8
29:00 252.5 F 16.7
30:00 252.4 F 16.6

Min 251.1 F 16.0
Max 252.7 F 16.8

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 177.1 F 0.0
5:00 199.5 F 0.0
10:00 201.3 F 0.0
15:00 208.1 F 0.0
20:00 218.5 F 0.0
25:00 226.3 F 0.0
30:00 233.0 F 0.0

Drying Complete

FILLING: 0:51
HEATING: 20:45
STERILIZING: 30:00
VENTING: 2:50
DRYING: 30:00
TOTAL CYCLE: 01:24:27
CYCLE COMPLETE

02/23/2021 10:45 AM
ACCEPTED: X
REJECTED:
Operator: U01