

-----  
Midmark M11 - v1.0.5  
Total Cycles: 2367

Name: RITTER2  
Sterilizer ID: V2313228  
05/15/2023 12:32 PM

BEGIN PACKS CYCLE  
Temp: 250 Degrees F  
Time: 30 Minutes  
Dry: 30 Minutes

FILLING CHAMBER

HEATING  
mm:ss Degrees PSI  
0:00 198.6 F 0.0  
2:00 194.9 F 4.1  
4:00 201.4 F 3.0  
6:00 216.7 F 5.6  
8:00 221.1 F 4.0  
10:00 226.3 F 5.2  
12:00 237.6 F 9.7  
14:00 250.1 F 15.9

STERILIZING  
mm:ss Degrees PSI  
0:00 251.1 F 16.6  
1:00 252.5 F 17.2  
2:00 252.2 F 16.9  
3:00 252.4 F 17.2  
4:00 251.9 F 17.3  
5:00 252.9 F 17.6  
6:00 252.8 F 17.4  
7:00 252.2 F 17.2  
8:00 251.7 F 17.0  
9:00 252.0 F 17.4  
10:00 253.1 F 17.7  
11:00 253.2 F 17.6  
12:00 253.0 F 17.6  
13:00 252.9 F 17.5  
14:00 252.7 F 17.4  
15:00 252.6 F 17.3  
16:00 252.4 F 17.3  
17:00 252.2 F 17.2  
18:00 251.9 F 17.1  
19:00 251.8 F 17.1  
20:00 251.7 F 17.0  
21:00 251.7 F 17.0  
22:00 251.7 F 17.1  
23:00 251.8 F 17.3  
24:00 251.9 F 17.4  
25:00 252.2 F 17.5  
26:00 252.5 F 17.7  
27:00 252.6 F 17.7  
28:00 253.3 F 17.7  
29:00 253.1 F 17.7  
30:00 253.0 F 17.7

Min 250.9 F 16.5  
Max 253.5 F 17.7

VENTING CHAMBER

DRYING  
Drying Start  
mm:ss Degrees PSI  
0:00 199.4 F 0.0  
5:00 199.5 F 0.0  
10:00 205.4 F 0.0  
15:00 205.8 F 0.0  
20:00 206.4 F 0.0  
25:00 206.4 F 0.0  
30:00 207.2 F 0.0  
Drying Complete

FILLING: 1:53  
HEATING: 14:54  
STERILIZING: 30:00

VENTING: 3:01  
DRYING: 30:00  
TOTAL CYCLE: 01:19:47  
CYCLE COMPLETE

05/15/2023 01:51 PM  
ACCEPTED:  
REJECTED:  
Reason: \_\_\_\_\_  
Operator: \_\_\_\_\_