
Midmark M11 - v1.0.0
Total Cycles: 864

Name: RITTER 1
Sterilizer ID: V2178198
06/07/2021 09:24 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 193.7 F 0.0
2:00 179.9 F 2.0
4:00 193.2 F 6.3
6:00 215.4 F 3.0
8:00 227.7 F 6.1
10:00 230.2 F 6.6
12:00 230.2 F 6.1
14:00 244.3 F 12.2

STERILIZING
mm:ss Degrees PSI
0:00 251.3 F 15.9
1:00 252.2 F 16.1
2:00 252.3 F 16.4
3:00 251.9 F 15.9
4:00 252.7 F 16.5
5:00 252.0 F 16.0
6:00 252.1 F 16.3
7:00 252.4 F 16.2
8:00 251.6 F 15.8
9:00 252.5 F 16.5
10:00 252.4 F 16.2
11:00 251.8 F 15.9
12:00 251.9 F 16.2
13:00 252.6 F 16.4
14:00 252.2 F 16.1
15:00 251.6 F 15.8
16:00 252.2 F 16.4
17:00 252.6 F 16.4
18:00 252.3 F 16.2
19:00 251.8 F 15.9
20:00 251.6 F 16.0
21:00 252.7 F 16.5
22:00 252.5 F 16.3
23:00 252.2 F 16.1
24:00 251.8 F 15.9
25:00 251.6 F 16.0
26:00 252.6 F 16.5
27:00 252.6 F 16.4
28:00 252.3 F 16.2
29:00 252.0 F 16.0
30:00 251.6 F 15.8

Min 251.3 F 15.8
Max 252.8 F 16.5

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 208.7 F 0.0
5:00 219.3 F 0.0
Dry Cycle Aborted

FILLING: 1:51
HEATING: 15:47
STERILIZING: 30:00
VENTING: 2:57
DRYING: 6:57
TOTAL CYCLE: 00:57:32

06/07/2021 10:21 AM

ACCEPTED: X
REJECTED:
Operator: U01