
Midmark M11 - v1.0.0
Total Cycles: 669

Name: MIDMARK
Sterilizer ID: V2178198
02/17/2021 08:44 AM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 188.2 F 0.0
2:00 179.5 F 1.6
4:00 180.4 F 5.1
6:00 205.8 F 1.9
8:00 216.9 F 5.0
10:00 222.3 F 3.9
12:00 232.4 F 7.4
14:00 231.3 F 6.4
16:00 243.9 F 11.9
18:00 251.0 F 15.9

STERILIZING
mm:ss Degrees PSI
0:00 251.6 F 16.1
1:00 252.1 F 16.2
2:00 252.6 F 16.5
3:00 251.4 F 16.0
4:00 252.1 F 16.4
5:00 251.3 F 16.1
6:00 252.3 F 16.5
7:00 251.5 F 16.1
8:00 252.2 F 16.6
9:00 252.3 F 16.3
10:00 251.4 F 16.0
11:00 252.5 F 16.6
12:00 252.0 F 16.3
13:00 251.4 F 16.0
14:00 252.2 F 16.6
15:00 252.1 F 16.5
16:00 251.8 F 16.2
17:00 251.3 F 16.1
18:00 252.4 F 16.6
19:00 252.2 F 16.5
20:00 251.8 F 16.2
21:00 251.2 F 16.0
22:00 252.0 F 16.6
23:00 252.0 F 16.5
24:00 251.9 F 16.3
25:00 251.4 F 16.1
26:00 251.7 F 16.4
27:00 252.0 F 16.6
28:00 252.0 F 16.4
29:00 251.5 F 16.2
30:00 251.2 F 16.1

Min 251.2 F 16.0
Max 252.8 F 16.6

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 173.3 F 0.0
5:00 187.4 F 0.0
10:00 184.2 F 0.0
15:00 204.5 F 0.0
20:00 194.8 F 0.0
25:00 164.3 F 0.0
Dry Cycle Aborted

FILLING: 1:18
HEATING: 18:37

STERILIZING: 30:00
VENTING: 3:05
DRYING: 27:26
TOTAL CYCLE: 01:20:25

02/17/2021 10:04 AM
ACCEPTED: X
REJECTED:
Operator: U01