

-----  
Midmark M11 - v1.0.0  
Total Cycles: 1025

Name: RITTER 1  
Sterilizer ID: V2178198  
07/19/2021 06:46 AM

BEGIN PACKS CYCLE  
Temp: 250 Degrees F  
Time: 30 Minutes  
Dry: 30 Minutes

FILLING CHAMBER

HEATING  
mm:ss Degrees PSI  
0:00 72.0 F 0.0  
2:00 81.4 F 0.0  
4:00 90.5 F 0.0  
6:00 85.6 F 0.0  
8:00 101.5 F 0.1  
10:00 130.6 F 0.1  
12:00 165.9 F 0.2  
14:00 196.9 F 1.0  
16:00 221.7 F 3.9  
18:00 228.0 F 6.4  
20:00 222.1 F 3.3  
22:00 233.1 F 7.3  
24:00 229.7 F 5.8  
26:00 226.6 F 4.7  
28:00 239.7 F 9.9  
30:00 250.3 F 15.0

STERILIZING  
mm:ss Degrees PSI  
0:00 251.3 F 15.7  
1:00 252.4 F 16.0  
2:00 251.5 F 16.0  
3:00 252.2 F 15.9  
4:00 251.4 F 16.0  
5:00 252.3 F 16.0  
6:00 251.2 F 15.6  
7:00 252.2 F 16.3  
8:00 252.3 F 16.0  
9:00 251.3 F 15.6  
10:00 251.8 F 16.2  
11:00 252.5 F 16.2  
12:00 252.2 F 15.9  
13:00 251.3 F 15.6  
14:00 251.5 F 16.1  
15:00 252.5 F 16.3  
16:00 252.5 F 16.1  
17:00 252.1 F 15.9  
18:00 251.4 F 15.7  
19:00 251.2 F 15.8  
20:00 252.1 F 16.3  
21:00 252.6 F 16.2  
22:00 252.5 F 16.1  
23:00 252.1 F 15.9  
24:00 251.5 F 15.7  
25:00 251.1 F 15.6  
26:00 251.3 F 16.0  
27:00 252.1 F 16.3  
28:00 252.6 F 16.3  
29:00 252.5 F 16.1  
30:00 252.4 F 16.0

Min 251.0 F 15.6  
Max 252.6 F 16.3

VENTING CHAMBER

DRYING  
Drying Start  
mm:ss Degrees PSI  
0:00 185.9 F 0.0  
5:00 194.8 F 0.0  
10:00 198.9 F 0.0  
15:00 198.1 F 0.0

20:00 198.2 F 0.0  
25:00 203.0 F 0.0  
30:00 173.2 F 0.0  
Drying Complete

FILLING: 0:50  
HEATING: 30:50  
STERILIZING: 30:00  
VENTING: 3:04  
DRYING: 30:00  
TOTAL CYCLE: 01:34:44  
CYCLE COMPLETE

07/19/2021 08:30 AM  
ACCEPTED: X  
REJECTED:  
Operator: U01