
Midmark M11 - v1.0.0
Total Cycles: 579

Name: MIDMARK
Sterilizer ID: V2178198
12/21/2020 02:41 PM

BEGIN PACKS CYCLE
Temp: 250 Degrees F
Time: 30 Minutes
Dry: 30 Minutes

FILLING CHAMBER

HEATING
mm:ss Degrees PSI
0:00 144.1 F 0.0
2:00 158.6 F 1.2
4:00 167.1 F 4.2
6:00 201.4 F 3.0
8:00 211.2 F 4.2
10:00 221.5 F 4.5
12:00 230.1 F 7.2
14:00 230.1 F 6.4
16:00 243.6 F 12.2

STERILIZING
mm:ss Degrees PSI
0:00 251.2 F 16.2
1:00 252.7 F 16.7
2:00 251.1 F 16.1
3:00 252.6 F 16.8
4:00 251.9 F 16.4
5:00 251.2 F 16.3
6:00 252.6 F 16.8
7:00 252.5 F 16.6
8:00 251.8 F 16.4
9:00 251.2 F 16.1
10:00 251.5 F 16.5
11:00 252.5 F 16.9
12:00 252.7 F 16.8
13:00 252.5 F 16.7
14:00 252.2 F 16.5
15:00 251.7 F 16.4
16:00 251.4 F 16.2
17:00 251.2 F 16.1
18:00 251.3 F 16.4
19:00 251.8 F 16.7
20:00 252.4 F 16.9
21:00 252.7 F 16.9
22:00 252.7 F 16.8
23:00 252.6 F 16.7
24:00 252.5 F 16.7
25:00 252.4 F 16.6
26:00 252.0 F 16.5
27:00 251.8 F 16.4
28:00 251.5 F 16.3
29:00 251.3 F 16.2
30:00 251.2 F 16.1

Min 251.1 F 16.1
Max 252.8 F 16.9

VENTING CHAMBER

DRYING
Drying Start
mm:ss Degrees PSI
0:00 193.7 F 0.0
5:00 203.2 F 0.0
10:00 207.3 F 0.0
15:00 209.5 F 0.0
20:00 208.3 F 0.0
25:00 212.4 F 0.0
30:00 204.4 F 0.0
Drying Complete

FILLING: 0:52
HEATING: 17:59

STERILIZING: 30:00
VENTING: 2:53
DRYING: 30:00
TOTAL CYCLE: 01:21:44
CYCLE COMPLETE

12/21/2020 04:02 PM
ACCEPTED: X
REJECTED:
Operator: U01